

Product Datasheet

TEX28 Antibody (1A8) H00001527-M05

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001527-M05

Updated 10/13/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001527-M05



H00001527-M05


TEX28 Antibody (1A8)

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1A8
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	1527
Gene Symbol	TEX28
Species	Human
Reactivity Notes	Human
Specificity/Sensitivity	TEX28 (1A8)
Immunogen	TEX28 (AAH26183, 1 a.a. ~ 411 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MVLKAEHTRSPSATLPSNVPSCRSLSSSEEDGPGSPSSLADGGGLAHLNLQD SVRHRILYLSEQLRVEKASRDGNTVSYLKLVSKADRHQVPHIQQAFQEKVNQRA SATIAQIEHRLHQCHQQLQELEEGCRPEGLLLMAESDPANCEPPSEKALLSEPP EPGGEDGPVNLPHASRPFILESRFQSLQGGTCLETEDVAQQQNLQKVKAELE EAKRFHISLQESYHSLKERSLTDLQLLLESLQEEKCRQALMEEQV
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	This antibody is useful for ELISA and Western Blot.




Images

Western Blot: TEX28 Antibody (1A8) [H00001527-M05] - Analysis of TEX28 expression in transfected 293T cell line by TEX28 monoclonal antibody (M05), clone 1A8. Lane 1: TEX28 transfected lysate (Predicted MW: 46.1 KDa). Lane 2: Non-transfected lysate.



Sandwich ELISA: TEX28 Antibody (1A8) [H00001527-M05] - Detection limit for recombinant GST tagged TEX28 is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00001527-M05

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NB600-986	Mouse IgG2a Kappa Light Chain Isotype Control (MG2a-53)
H00001527-P01-10ug	TEX28 Recombinant Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001527-M05

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

