

Product Datasheet

TEX28 Antibody (1A8)

H00001527-M05

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001527-M05

Updated 12/26/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001527-M05



H00001527-M05


TEX28 Antibody (1A8)

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1A8
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	1527
Gene Symbol	TEX28
Species	Human
Specificity/Sensitivity	TEX28 (1A8)
Immunogen	TEX28 (AAH26183, 1 a.a. - 410 a.a.) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MVLKAEHTRSPSATLPSNVPSCRSLSSSEEDGPSGSSLADGGGLAHNLQDSVR HRILYLSEQLRVEKASRDGNTVSYLKLVSADRHRQVPHIQQAFKVNQRASATI AQIEHRLHQCHQQQLQELEEGCRPEGLLLMAESDPANCEPPSEKALLSEPPEPG GEDGPVNLPHASRPFILESRFQSLQQGTCLLETEDVAQQQNLLLQKVKAEELEEA KRFHISLQESYHSLKERSLTDLQLLLESLQEEKCRQALMEEQVNGRLQGQLNEI YNLKHNLACSEERMAYLSYERAKEIWEITETFKSRISKLEMLQQVTQLEAAEHL QSRPPQMLFKFLSPRLSLATVLLVFVSTLCACPSSLISSRLCTCTMLMLIGLV AWQRWRAIPATDWQEWVPSRCRLYSKDSGPPADGP
Notes	Quality control test: Antibody Reactive Against Recombinant Protein. This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot, ELISA
Application Notes	It has been used for ELISA and WB.




Images

Western Blot: TEX28 Antibody (1A8) [H00001527-M05] - Analysis of TEX28 expression in transfected 293T cell line by TEX28 monoclonal antibody (M05), clone 1A8. Lane 1: TEX28 transfected lysate (Predicted MW: 46.1 KDa). Lane 2: Non-transfected lysate.



Sandwich ELISA: TEX28 Antibody (1A8) [H00001527-M05] - Detection limit for recombinant GST tagged TEX28 is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00001527-M05

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43318-0.5mg	Mouse IgG2a Kappa Light Chain Isotype Control (M2a)
H00001527-P01-10ug	Recombinant Human TEX28 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001527-M05

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

