

Product Datasheet

DAB2 Antibody H00001601-B01P

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001601-B01P

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001601-B01P



H00001601-B01P

DAB2 Antibody

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Mouse
Gene ID	1601
Gene Symbol	DAB2
Species	Human
Immunogen	DAB2 (AAH03064.1, 1 a.a. - 770 a.a.) full-length human protein. MSNEVETSATNGQPQQAAAPKAPSKKEKKKGPEKTDEYLLARFKGDGVKYKA KLGIDDDVPDARGDKMSQDSMMKLKGMAGRSQGGQHKQRIWVNIISLGGIKIID EKTGVIEHEHPVNIKISFIARDVTDNRAFGYVCGGEGQHQQFFAIKTGQQAEPLVV DLKDLFQVIYNVKKKEEKKKIEEASKAVENGSEALMILDDQTNKLKSGVDQMD LFGDMSTPPDLNSPTESKDILLVDLNSEIDTNQNSLRENPFNTGITSCLPRPT PQASFLPENAFSANLNFFPTPNPDPFRDDPFTQPDQSTPSSFDSLKSPDQKKE NSSSSSTPLSNGPLNGDVDYFGQQFDQISNRTGKQEAQAGPWPFFSSQTQPA VRTQNGVSEREQNGFSVKSSPNPFVGSPPKGLSIQNGVKQDLESSVQSSPHD SIAIIPPPQSTKPRGRRTAKSSANDLLASDIFAPPVSEPSGQASPTGQPTALQP NPLDLFKTSAPAPVGPLVGLGGVTVTLPPQAGPWNTASLVFNQSPSMAPGAMM GGQPSGFSQPVIFGTSPAVSGWNQPSPFPAASTPPPVPVWVWGPSASVAPNAW STTSPLGNPFQSNIFPAPAVSTQPPSMHSSLLVTPPQPPPRAGPPKDISSDAFT ALDPLGDKEIKDVKEMFKDFQLRQPPAVPARKGEQTSSGTLFAFASYFNSKVGI PQENADHDDFDANQLLNKINEPPKPAPRQVSLPVTKSTDNAFENPFFKDSFGS SQASVASSQPVSSEMYRDPFGNPFA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/Immunofluorescence
Application Notes	It has been used for WB and IF.




Images


Western Blot: DAB2 Antibody [H00001601-B01P] - Analysis of DAB2 expression in human liver.




Immunocytochemistry/Immunofluorescence: DAB2 Antibody [H00001601-B01P] - Analysis of purified antibody to DAB2 on HeLa cell. (antibody concentration 10 ug/ml)



Western Blot: DAB2 Antibody [H00001601-B01P] - Analysis of DAB2 expression in HeLa.



Western Blot: DAB2 Antibody [H00001601-B01P] - Analysis of DAB2 expression in transfected 293T cell line by DAB2 polyclonal antibody. Lane 1: DAB2 transfected lysate(84.7 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001601-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

