

Product Datasheet

Desmocollin-2 Antibody H00001824-B01P

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001824-B01P

Updated 8/29/2017 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001824-B01P



H00001824-B01P

Desmocollin-2 Antibody


Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)

Product Description	
Host	Mouse
Gene ID	1824
Gene Symbol	DSC2
Species	Human
Immunogen	DSC2 (NP_004940.1, 1 a.a. - 847 a.a.) full-length human protein. MEAARPSGSWNGALCRLLLTLAILIFASDACKNVTLHVPSKLDKAEKLVGRVNL KECFTAANLIHSSDPDFQILEDGSVYTTNTILLSSEKRSFTILLSNTENQEKKKIFV FLEHQTKVLKKRHTKEKVLRRAKRRWAPIPCSMLENSLGPFLFLQQVQSDTA QNYTIYYSIRGPGVDQEPRNLFYVERDTGNLYCTRVPDREQYESFEIIFATTPD GYTPELPLPLIKIEDENDNYPIFTEETYFTIFENCVRGTTVGQVCATDKDEPDT MHTRLKYSIIGQVPPSPTLFSMHPTTGVIITSSQLDRELIDKYQLKIKVQDMDG QYFGLQTTSTCIINIDDVNDHLPFTFRTSYVTSVEENTVDVEILRVTVEDKDLVNT ANWRANYTILKGNENGNFKIVTDAKTNEGVLCVVKPLNYEEKQQMILQIGVVNE APFSREASPRSAMSTATVTNVVEDQDEGPECNPPIQTVRMKENAEVGTTSNG YKAYDPETRSSSGIRYKCLTDPTGWVTIDENTGSIKVFRLDREAETIKNGIYNIT VLASDQGGRTCTGTLGILQDVNDNSPFIPKKTVIICKPTMSSAEIVAVDPDEPIH GPPFDLSLESSTSEVQRMWRLKAINDTAARLSYQNDPPFGSYVVPITVRDRLG MSSVTSLDVTLCDCITENDCTHRVDPRIGGGGVQLGKWAILLAILLFCILFT LVCGASGTSKQPKVIPDDLAQQNLIVSNTEAPGDDKVYSANGFTTQTVGASAQ GVCQTVGSGIKNGGQETIEMVKGGHQTSESCRAGHHHTLDSCRGGHTEVD NCRYTYSEWHSFTQPRLGEESIRGHTLIKN
Notes	Quality control test: Antibody reactive against mammalian transfected lysate. This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	It has been used for WB.

Images

Western Blot: Desmocollin-2 Antibody [H00001824-B01P] - Analysis of DSC2 expression in transfected 293T cell line by DSC2 polyclonal antibody. Lane 1: DSC2 transfected lysate(93.17 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00001824-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBP1-85710PEP	Desmocollin-2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001824-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

