

# Product Datasheet

## DUSP9 Antibody - Azide and BSA Free H00001852-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00001852-D01P](http://www.novusbio.com/H00001852-D01P)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00001852-D01P](http://www.novusbio.com/reviews/destination/H00001852-D01P)



**H00001852-D01P**

DUSP9 Antibody - Azide and BSA Free

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 mg                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | No Preservative                                                                                                  |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Protein G purified                                                                                               |
| <b>Buffer</b>        | PBS (pH 7.4)                                                                                                     |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Quality control test: Antibody reactive against mammalian transfected lysate.                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Gene ID</b>     | 1852                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Gene Symbol</b> | DUSP9                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Species</b>     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen</b>   | DUSP9 (AAH60837.1, 1 a.a. - 384 a.a.) full-length human protein.<br>MEGLGRSCLWLRRELSPPRPRLLLDCRSRELYESARIGGALSVALPALLRRL<br>RRGSLSVRALLPGPPLQPPPPAPVLLYDQGGGRRRRGEAEAEAEWEAESVL<br>GTLLQKLREEGYLAYYLQGGFSRFQAECPHLCETSLAGRAGSSMAPLPGPVP<br>VVGLGSLCLGSDCSDAESEADRDSMSCGLDSEGATPPPVGRLASFVPQILPNL<br>YLGSARDSANLESLAKLGIRYILNVTNLPNFFFEKNGDFHYKQIPISDHWSQNLS<br>RFFPEAIEFIDEALSQNRGVLVHCLAGVSRSVTVTVAYLMQKLHLSLNDAYDLV<br>KRKKSNI SPNFNFMGQLLDFERSLRLEERHSQE QGSGGQASAASNPPSFFFTTP<br>TSDGAFELAPT |
| <b>Notes</b>       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                                                                                               |

**Product Application Details**

|                              |                          |
|------------------------------|--------------------------|
| <b>Applications</b>          | Western Blot             |
| <b>Recommended Dilutions</b> | Western Blot             |
| <b>Application Notes</b>     | It has been used for WB. |

**Images**

Western Blot: DUSP9 Antibody [H00001852-D01P] - Analysis of DUSP9 expression in transfected 293T cell line by DUSP9 polyclonal antibody. Lane 1: DUSP9 transfected lysate (41.90 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00001852-D01P**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00001852-P01-10ug | Recombinant Human DUSP9 GST (N-Term) Protein        |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00001852-D01P](http://www.novusbio.com/reviews/submit/H00001852-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

