

Product Datasheet

Ephrin-A3 Antibody (2H3) - Azide and BSA Free H00001944-M10

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00001944-M10

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00001944-M10



H00001944-M10

Ephrin-A3 Antibody (2H3) - Azide and BSA Free

Product Information

Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2H3
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description


Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	1944
Gene Symbol	EFNA3
Species	Human
Specificity/Sensitivity	EFNA3 - ephrin-A3 (2H3)
Immunogen	EFNA3 (NP_004943, 45 a.a. ~ 145 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. RREGYTVQVNVNDYLDIYCPHYNSSGVGPGAGPGPGGGAEQYVLYMVS RNG YRTCNASQGFKRWE CNRPHAPHSPIKFSEKFQRYSAFSLGYEFHAGHEY Y
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details

Applications	ELISA, Immunoprecipitation
Recommended Dilutions	ELISA, Immunoprecipitation
Application Notes	It has been used for ELISA.

Images

Immunoprecipitation: Ephrin-A3 Antibody (2H3) [H00001944-M10] - Analysis of EFNA3 transfected lysate using anti-EFNA3 monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with EFNA3 rabbit polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00001944-M10

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP2-77495-50ug	Recombinant Human Ephrin-A3 hIgG-His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00001944-M10

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

