

# Product Datasheet

## Flt-3 Ligand/FLT3L Antibody

### H00002323-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00002323-D01P](http://www.novusbio.com/H00002323-D01P)

Updated 6/27/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00002323-D01P](http://www.novusbio.com/reviews/destination/H00002323-D01P)



**H00002323-D01P**

Flt-3 Ligand/FLT3L Antibody


| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 mg                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | No Preservative                                                                                                  |
| Isotype             | IgG                                                                                                              |
| Purity              | Protein G purified                                                                                               |
| Buffer              | PBS (pH 7.4)                                                                                                     |

| Product Description |                                                                                                                                                                                                                                                                                                                                  |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Quality control test: Antibody reactive against mammalian transfected lysate.                                                                                                                                                                                                                                                    |
| Host                | Rabbit                                                                                                                                                                                                                                                                                                                           |
| Gene ID             | 2323                                                                                                                                                                                                                                                                                                                             |
| Gene Symbol         | FLT3LG                                                                                                                                                                                                                                                                                                                           |
| Species             | Human                                                                                                                                                                                                                                                                                                                            |
| Immunogen           | FLT3LG (NP_001450.2, 1 a.a. - 235 a.a.) full-length human protein.<br>MTVLAPAWSPPTYLLLLLLLLSSGLSGTQDCSFQHSPISDFAVKIRELSDYLLQD<br>YPVTVASNLQDEELCGGLWRLVLAQRWMERLKT VAGSKMQGLLERVNTIHF<br>VTKCAFQPPPSCLRFVQTNISRLQETSEQLVALKPWITRQNF SRCLELQCQPD<br>SSTLPPPWSRPLEATAPTAPQPPLLLLLLLPVGLLLLAAAWCLHWQRTRRRTP<br>RPGEQVPPVPSPQDLLLVEH |
| Notes               | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                               |

| Product Application Details |                          |
|-----------------------------|--------------------------|
| Applications                | Western Blot             |
| Recommended Dilutions       | Western Blot             |
| Application Notes           | It has been used for WB. |

**Images**

Western Blot: Flt-3 Ligand Antibody [H00002323-D01P] - Analysis of FLT3LG expression in transfected 293T cell line by FLT3LG polyclonal antibody. Lane 1: FLT3LG transfected lysate (26.40 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis.  
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00002323-D01P](http://www.novusbio.com/reviews/submit/H00002323-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

