

Product Datasheet

HHEX Antibody (1A5) - Azide and BSA Free H00003087-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00003087-M02

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00003087-M02



H00003087-M02


HHEX Antibody (1A5) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1A5
Preservative	No Preservative
Isotype	IgG2a Lambda
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	3087
Gene Symbol	HHEX
Species	Human
Specificity/Sensitivity	HHEX - homeobox, hematopoietically expressed
Immunogen	HHEX (AAH14336, 1 a.a. ~ 270 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MQYPHPGPAAGAVGVPLYAPTPLLQPAHPTPFYIEDILGRGPAAPTPAPTLPSP NSSFTSLVSPYRTPVYEPTPIHPAFSHHSAALAAAYGPGGFGGPLYPFPRTVN DYTHALLRHDPLGKPLLWSPFLQRPLHKRKGQVRFSDNDQTIELEKKFETQKY LSPPERKRLTKMLQLSERQVKTWFQNRRAKWRLKQENPQSNKKEELESLS SCDQRQDLPSEQNKGASLDSSQCSPSPASQEDLESEISEDSDQEVDIEGDKSY FNAG
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA
Recommended Dilutions	Western Blot 1:500, ELISA, Immunocytochemistry/ Immunofluorescence, Sandwich ELISA
Application Notes	Antibody reactivity against cell lysate, tissue lysate, transfected lysate and recombinant protein for WB. It has been used for IF and ELISA.




Images


Western Blot: HHEX Antibody (1A5) [H00003087-M02] - Analysis of HHEX expression in transfected 293T cell line by HHEX monoclonal antibody (M02), clone 1A5. Lane 1: HHEX transfected lysate (Predicted MW: 30 KDa). Lane 2: Non-transfected lysate.




Immunocytochemistry/Immunofluorescence: HHEX Antibody (1A5) [H00003087-M02] - Immunofluorescence of monoclonal antibody to HHEX (H00003087-M02) on HeLa cell . [antibody concentration 10 ug/ml]




Western Blot: HHEX Antibody (1A5) [H00003087-M02] - Western Blot analysis of HHEX expression in MES-SA/Dx5 (Cat # L021V1).




Western Blot: HHEX Antibody (1A5) [H00003087-M02] - Western Blot analysis of HHEX expression in SW-13 (Cat # L005V1).



Western Blot: HHEX Antibody (1A5) [H00003087-M02] - Western Blot analysis of HHEX expression in Hela NE (Cat # L013V3).




Western Blot: HHEX Antibody (1A5) [H00003087-M02] - Western Blot analysis of HHEX expression in H9c2(2-1).




Western Blot: HHEX Antibody (1A5) [H00003087-M02] - Western Blot analysis of HHEX expression in human liver.



Western Blot: HHEX Antibody (1A5) [H00003087-M02] - HHEX monoclonal antibody (M02), clone 1A5. Analysis of HHEX expression in human placenta.



Sandwich ELISA: HHEX Antibody (1A5) [H00003087-M02] - Detection limit for recombinant GST tagged HHEX is approximately 0.3ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00003087-M02

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP2-23108	Recombinant Human HHEX His Protein
2496-ZN-010	Insulysin/IDE [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00003087-M02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

