

Product Datasheet

IL-13R alpha 1 Antibody H00003597-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00003597-B01P

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00003597-B01P



H00003597-B01P**IL-13R alpha 1 Antibody**

Product Information	
Unit Size	0.05 mg
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Unpurified
Buffer	PBS (pH 7.4)


Product Description	
Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Mouse
Gene ID	3597
Gene Symbol	IL13RA1
Species	Human
Specificity/Sensitivity	IL13RA1 - interleukin 13 receptor, alpha 1,
Immunogen	IL13RA1 (NP_001551.1, 1 a.a. - 427 a.a.) full-length human protein. MEWPARLCGLWALLLCAGGGGGGGGAAPTETQPPVTNLSVSVENLCTVIWT WNPPEGASSNCSLWYFSHFQDKQDKKIAPETRRSIEVPLNERICLQVGSQCST NESEKPSILVEKCISPPEGDPESAVTELQCIWHNLSYMKCSWLPGRNTSPDTN YTLYYWHRSLKIHQCENIFREGQYFGCSFDLTQVKDSSFEQHSVQIMVKDNA GKIKPSFNIVPLTSRVKPDPPHIKNLSFHNDLDYVQWENPQNFISRCLFYEEVN NSQTETHNVFYVQEAKCENPEFERNVENTSCFMVPGVLPDTLNTVRIRVKTNK LCYEDDKLWSNWSQEMSIGKKRNSTLYITMLLIVPVIVAGAIIVLLLYLKRLKIIIFP PIPDGPKIFKEMFGDQNDLTLHWKKYDIYEKQTKKEETDSVVLIENTLKKASQ
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Flow Cytometry
Recommended Dilutions	Western Blot, Flow Cytometry
Application Notes	This antibody is reactive against transfected lysate in WB, FACS analysis and as a detection antibody in ELISA. It has been use for FACS.




Images

Western Blot: IL-13 R alpha 1 Antibody [H00003597-B01P] - Analysis of IL13RA1 expression in transfected 293T cell line by IL13RA1 polyclonal antibody. Lane 1: IL13RA1 transfected lysate(46.97 KDa). Lane 2: Non-transfected lysate.



Flow Cytometry: IL-13 R alpha 1 Antibody [H00003597-B01P] - Analysis of negative control 293 cells (Black) and IL13RA1 expressing 293 cells (Green) using IL13RA1 purified mouse polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00003597-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

