

Product Datasheet

Kininogen Antibody H00003827-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00003827-B01P

Updated 7/16/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00003827-B01P



H00003827-B01P

Kininogen Antibody

| Product Information | |
|---------------------|--|
| Unit Size | 0.05 mg |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | No Preservative |
| Isotype | IgG |
| Purity | Protein A purified |
| Buffer | PBS (pH 7.4) |

| Product Description | |
|---------------------|--|
| Description | Quality control test: Antibody reactive against mammalian transfected lysate. |
| Host | Mouse |
| Gene ID | 3827 |
| Gene Symbol | KNG1 |
| Species | Human |
| Immunogen | Kininogen (NP_000884.1, 1 a.a. - 427 a.a.) full-length human protein. MKLITILFLCSRLLLLSLTQESQSEEIDCNDKDLFKAVDAALKKYNSQNSQNSNNQFV LYRITEATKTVGSDTFYSFKYEIKEGDCPVQSGKTWQDCEYKDAKAATGECT ATVGKRSSTKFSVATQTCQITPAEGPVVTAQYDCLGCVHPISTQSPDLEPILRH GIQYFNNNTQHSSLFMLNEVKRAQRQVVAGLNFRITYSIVQTNCSENFLFLTP DCKSLWNGDTGECTDNAYIDIQLRIASFQNCDIYPGKDFVQPPTKICVGCPRDI PTNSPELEETLHTITKLNAENNATFYFKIDNVKKARVQVVAGKKYFIDFVARET TCSKESNEELTESCETKKLQSLDCNAEVVVPWEKKIYPTVNCQPLGMISLM KRPPGFSPFRSSRIGEIKEETTSHLRSCYKGRPPKAGAEPASEREVS |
| Notes | This product is produced by and distributed for Abnova, a company based in Taiwan. |


| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | Antibody reactivity against Recombinant Protein with GST tag on ELISA and WB and also on transfected lysate in WB. GST tag alone is used as a negative control. |

Images

Western Blot: Kininogen Antibody [H00003827-B01P] - Analysis of Kininogen expression in human pancreas.



Western Blot: Kininogen Antibody [H00003827-B01P] - Analysis of Kininogen expression in transfected 293T cell line by Kininogen polyclonal antibody. Lane 1: Kininogen transfected lysate(46.97 KDa). Lane 2: Non-transfected lysate.



Publications

Abdullah-Soheimi SS, Lim BK, Hashim OH et al. Patients with ovarian carcinoma excrete different altered levels of urine CD59, kininogen-1 and fragments of inter-alpha-trypsin inhibitor heavy chain H4 and albumin. *Proteome Sci*;8:58. 2010-11-17 [PMID: 21083881]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00003827-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

