

Product Datasheet

MSX2 Antibody H00004488-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00004488-B01P

Updated 8/29/2017 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00004488-B01P



H00004488-B01P

MSX2 Antibody

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	4488
Gene Symbol	MSX2
Species	Human
Immunogen	MSX2 (AAH15509.1, 1 a.a. - 267 a.a.) full-length human protein. MASPSKGNLDFSPDEEGPAVVAGPGPGGGAEGAAEERRVKVSSLPFSVEAL MSDKKPPKEASPLPAESASAGATLRPLLLSGHGAREAHSPGPLVKPFETASVK SENSEDGAAWMQEPGRYSPPRHTSPTTCTLRKHKTNRKPRTPTTSSQLLAL ERKFRQKQYLSIAERAEFSSSLNLTETQVKIWFQNRRAKAKRLQEAELEKLKMA AKPMLPSSFSLPFISSPLQAASIYGASYPFHRPVLPIPPVGLYATPVGYGMYHL S
Notes	Quality control test: Antibody reactive against mammalian transfected lysate. This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/Immunofluorescence
Application Notes	It has been used for WB and Functional.




Images

Western Blot: MSX2 Antibody [H00004488-B01P] - Analysis of MSX2 expression in transfected 293T cell line by MSX2 polyclonal antibody. Lane 1: MSX2 transfected lysate(29.37 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: MSX2 Antibody [H00004488-B01P] - Analysis of purified antibody to MSX2 on HeLa cell. (antibody concentration 10 ug/ml)





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00004488-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00004488-Q03-10ug	Recombinant Human MSX2 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00004488-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

