

Product Datasheet

PCM1 Antibody H00005108-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00005108-B01P

Updated 8/29/2017 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00005108-B01P




H00005108-B01P

PCM1 Antibody

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	5108
Gene Symbol	PCM1
Species	Human
Specificity/Sensitivity	PCM1 - pericentriolar material 1,
Immunogen	PCM1 (AAH00453.1, 1 a.a. - 530 a.a.) full-length human protein. MATGGGPFEDGMNDQDLPNWSNENVDDRLNNMDWGAQQKKANRSSEKNKK KFGVESDKRVTNDISPESPGVGRRRTKTPHTFPHSRYSQMSVPEQAELEK LKQRINFSDLDQRSIGSDSQGRATAANNKRQLSENRPFNFLPMQINTNKS ASTSPPNRETIGSAQCKELFASALSNDLLQNCQVSEEDGRGEPAMESSQIVSR LVQIRDYITKASSMREDLVEKNERSANVERLTHLIDHLKEQEKSVMKFLKKILAR ENEEEDVRTIDSAVGSGSVAESTSLNIDVQSEASDTTARDPQQEPMEEIENLKK QHDLLKRMLQQQEQLRALQGRQAALLALQHKAEQAIAMDDSVVAETAGSLS GVSITSELNEELNDLIQRFHNQLRDSQPPAVPDNRRQAESLSLTREVSQSRKP SASERLPDEKVELFSKMRVLQEKKQKMDKLLGELHTLRDQHLNNSSSSPQRS VDQRSTSAPSASVGLAPVVNGESNSLTSSVPYPTASLVSQNESENEGHLNPSE KLQ
Notes	Quality control test: Antibody reactive against mammalian transfected lysate. This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Antibody reactivity against Recombinant Protein with GST tag on WB and also on transfected lysate in WB. GST tag alone is used as a negative control.

Images

Western Blot: PCM1 Antibody [H00005108-B01P] - Analysis of PCM1 expression in transfected 293T cell line by PCM1 polyclonal antibody. Lane 1: PCM1 transfected lysate(58.3 KDa). Lane 2: Non-transfected lysate.



Publications

Akimov V, Rigbolt KT, Nielsen MM, Blagoev B. Characterization of ubiquitination dependent dynamics in growth factor receptor signaling by quantitative proteomics. *Mol Biosyst.* 2011 Sep 28 [PMID: 21956701]



Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00005108-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00005108-P01-10ug	Recombinant Human PCM1 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00005108-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

