

Product Datasheet

PPP2R2B Antibody (1B10) - Azide and BSA Free H00005521-M23

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00005521-M23

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00005521-M23



H00005521-M23


PPP2R2B Antibody (1B10) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1B10
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	5521
Gene Symbol	PPP2R2B
Species	Human
Specificity/Sensitivity	PPP2R2B - protein phosphatase 2 (formerly 2A), regulatory subunit B, beta isoform (1B10)
Immunogen	PPP2R2B (AAH31790, 1 a.a. ~ 443 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MEEDIDTRKINNSFLRDHSYATEADIISTVEFNHTGELLATGDKGGRVVIFQREQ ESKNQVHRRGEYNVYSTFQSHEPEFDYLSLEIEEKINKIRWLPQQNAAYFLLS TNDKTVKLWKVSEDRKPEGYNLKDEEGRLRDPATITTLRVPVLRPMDLMVEA TPRRVFANAHTYHINSISVNSDYETYMSADDLRINLWNFEITNQSFNIVDIKPAN MEELTEVITAAEFHPHCNTFVYSSSKGTIRLCDMRASALCDRHTKFFEEPEDP SNRSFFSEIISISDVKFSGRYIMTRDYLTVKVWDLNMENRPIETYQVHDYLR SKLCSLYENDCIFDKFECVWNGSDSVIMTGSYNFFRMFDRNTKRDVTLEASR ENSKPRAILKPRKVCVGGKRRKDEISVDSLDFSKKILHTAWHPSENIIVAATNN LYIFQDKVN
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	ELISA, Immunoprecipitation
Recommended Dilutions	ELISA, Immunoprecipitation
Application Notes	It has been used for ELISA.



Images

Immunoprecipitation: PPP2R2B Antibody (1B10) [H00005521-M23] - Analysis of PPP2R2B transfected lysate using anti-PPP2R2B monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with PPP2R2B rabbit polyclonal antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00005521-M23

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
H00005521-P01-10ug	Recombinant Human PPP2R2B GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00005521-M23

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

