Product Datasheet

PTMS Antibody (3H5) - Azide and BSA Free H00005763-M17-100ug

Unit Size: 100 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/H00005763-M17-100ug

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/H00005763-M17-100ug



H00005763-M17-100ug


PTMS Antibody (3H5) - Azide and BSA Free	
Product Information	
Unit Size	100 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3H5
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse PTMS Antibody (3H5) - Azide and BSA Free (H00005763-M17-100ug) is a monoclonal antibody validated for use in WB,
	ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	J ,
Host Gene ID	guarantee.
	guarantee. Mouse
Gene ID	guarantee. Mouse 5763
Gene ID Gene Symbol	guarantee. Mouse 5763 PTMS
Gene ID Gene Symbol Species	guarantee. Mouse 5763 PTMS Human PTMS (NP_002815.3, 1 a.a. ~ 102 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MSEKSVEAAAELSAKDLKEKKEKVEEKASRKERKKEVVEEEENGAEEEEETA

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence







Western Blot: PTMS Antibody (3H5) [H00005763-M17-100ug] - Detection against Immunogen (37.9 KDa).




Immunocytochemistry/Immunofluorescence: PTMS Antibody (3H5) [H00005763-M17-100ug] - Monoclonal antibody to PTMS on HeLa cell . [antibody concentration 10 ug/ml]



Western Blot: PTMS Antibody (3H5) [H00005763-M17-100ug] - Analysis of PTMS expression in transfected 293T cell line by PTMS monoclonal antibody (M17), clone 3H5. Lane 1: PTMS transfected lysate (Predicted MW: 11.5 KDa). Lane 2: Non-transfected lysate.



ELISA: PTMS Antibody (3H5) [H00005763-M17-100ug] - Detection limit for recombinant GST tagged PTMS is 0.3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112

USA

Phone: 303.730.1950 Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6

Canada

Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-

techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to H00005763-M17-100ug

HAF007 Goat anti-Mouse IgG Secondary Antibody [HRP]

NB720-B Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]

NBP1-96981-0.5mg Mouse IgG2a Kappa Isotype Control (M2AK)

NBP2-23316 Recombinant Human PTMS His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00005763-M17-100ug

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

