

# Product Datasheet

## SSB Antibody H00006741-B01P-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00006741-B01P-50ug](http://www.novusbio.com/H00006741-B01P-50ug)

Updated 1/12/2021 v.20.1

### Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00006741-B01P-50ug](http://www.novusbio.com/reviews/destination/H00006741-B01P-50ug)



**H00006741-B01P-50ug**

SSB Antibody

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	6741
Gene Symbol	SSB
Species	Human
Immunogen	SSB (ABM85385.1, 1 a.a. - 408 a.a.) full-length human protein. MAENGDNEKMAALEAKICHQIEYYFGDFNLPRDKFLKEQIKLDEGWVPLEIMIKF NRLNRLTTDFNVIVEALSLSKAELMEISEDKTKIRRSPSKPLPEVTDEYKNDVKN RSVYIKGFPTDATLDDIKEWLEDKGQVLNIQMRRTLHKAFKGSIFVVFDSIESAK KFVETPGQKYKETDLLILFKDDYFAKKNEERKQNKVEAKLRAKQEQAQKLE EDAEMKSLEEKIGCLLKFSGDLDDQTCREDLHILFSNHGEIKWIDFVRGAKEGIL FKEKAKEALGKAKDANNGNLQLRNKEVTWEVLEGEVEKEALKKIIEDQQESLN KWKSKGRRFKGKGKGNKAAQPGSGKGKVQFQGKTKFASDDEHDEHDENG ATGPVKRAREETDKEEPASKQQKTENGAGDQ
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/Immunofluorescence




## Images

Western Blot: SSB Antibody [H00006741-B01P-50ug] - Analysis of SSB expression in transfected 293T cell line (H00006741-T01) by SSB MaxPab polyclonal antibody. Lane 1: SSB transfected lysate(44.88 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: SSB Antibody [H00006741-B01P-50ug] - Purified MaxPab antibody to SSB on HeLa cell. [antibody concentration 10 ug/ml]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00006741-B01P-50ug](http://www.novusbio.com/reviews/submit/H00006741-B01P-50ug)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

