

Product Datasheet

Elongin A Antibody (1F3) - Azide and BSA Free H00006924-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00006924-M02

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00006924-M02



H00006924-M02

Elongin A Antibody (1F3) - Azide and BSA Free

Product Information

Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1F3
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4

Product Description

Description	Quality control test: Antibody Reactive Against Recombinant Protein.
Host	Mouse
Gene ID	6924
Gene Symbol	TCEB3
Species	Human
Specificity/Sensitivity	TCEB3 (1F3)
Immunogen	TCEB3 (NP_003189, 81 a.a. ~ 190 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. NAEPDEQDFEKSNSRKRPRDALQKEEEMEGDYQETWKATGSRSPDHRQK KHRKLSELERPHKVSHGHERDERKRCHRMSPPTYSSDPESSDYGHVQSPPS CTSPHQMY
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details

Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Application Notes	Antibody reactivity against cell lysate and recombinant protein for WB. It has also been used for IF and ELISA.




Images


Western Blot: Elongin A Antibody (1F3) [H00006924-M02] - TCEB3 monoclonal antibody (M02), clone 1F3 Analysis of TCEB3 expression in Hela S3 NE.



Immunocytochemistry/Immunofluorescence: Elongin A Antibody (1F3) [H00006924-M02] - Analysis of monoclonal antibody to TCEB3 on HeLa cell. Antibody concentration 10 ug/ml.



Western Blot: Elongin A Antibody (1F3) [H00006924-M02] - Analysis of TCEB3 expression in transfected 293T cell line by TCEB3 monoclonal antibody (M02), clone 1F3. Lane 1: TCEB3 transfected lysate(87.2 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00006924-M02

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP1-87040PEP	Elongin A Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00006924-M02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

