

# Product Datasheet

## Tensin 1 293T Cell Transient Overexpression Lysate H00007145-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00007145-T01](http://www.novusbio.com/H00007145-T01)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00007145-T01](http://www.novusbio.com/reviews/destination/H00007145-T01)



**H00007145-T01****Tensin 1 293T Cell Transient Overexpression Lysate****Product Information**

|                      |   |
|----------------------|---|
| <b>Unit Size</b>     | 0.1 ml  |
| <b>Concentration</b> | mg/ml   |
| <b>Storage</b>       | Store at -80C. Avoid freeze-thaw cycles.  |
| <b>Buffer</b>        | 1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue) |

**Product Description**


|                                |   |
|--------------------------------|---|
| <b>Description</b>             | Quality control test: Transient overexpression cell lysate was tested with Anti-TNS1 antibody by Western Blots. <b>Plasmid:</b> pCMV-TNS1 full-length |
| <b>Gene ID</b>                 | 7145  |
| <b>Gene Symbol</b>             | TNS1  |
| <b>Species</b>                 | Human   |
| <b>Specificity/Sensitivity</b> | Tensin 1 293T Cell Transient Overexpression Lysate(Denatured)   |
| <b>Notes</b>                   | This product is produced by and distributed for Abnova, a company based in Taiwan.  |
| <b>Lysate Type</b>             | 293T Cell Transient Overexpression  |
| <b>Lysate Protein State</b>    | Denatured   |

**Product Application Details**

|                              |              |
|------------------------------|--------------|
| <b>Applications</b>          | Western Blot |
| <b>Recommended Dilutions</b> | Western Blot |

**Images**

Lane 1: TNS1 transfected lysate (43.7 KDa) Lane 2: Non-transfected lysate.



TNS1 transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00007145-T01**

---

|                    |  |
|--------------------|--|
| H00007145-P01-10ug | Recombinant Human Tensin 1 GST (N-Term) Protein  |
| 210-TA-005         | TNF-alpha [Unconjugated]                         |
| NBP1-84129         | Tensin 1 Antibody                                |
| MAB1455            | Albumin Antibody (188835) [Unconjugated] - Serum |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00007145-T01](http://www.novusbio.com/reviews/submit/H00007145-T01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

