

Product Datasheet

Frizzled-5 Antibody H00007855-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00007855-D01P

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00007855-D01P



H00007855-D01P**Frizzled-5 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Quality control test: Antibody reactive against mammalian transfected lysate.
Host	Rabbit
Gene ID	7855
Gene Symbol	FZD5
Species	Human, Mouse
Specificity/Sensitivity	Reacts with frizzled homolog 5 (Drosophila).
Immunogen	FZD5 (ENSP00000354607, 1 a.a. - 585 a.a.) full-length human protein. MARPDPSAPPSLLLLLLAQLVGRAAAASKAPVCQEITVPMCRGIGYNLTHMPN QFNHDTQDEAGLEVHQFWPLVEIQCSDDLRFLLCSMYTPICLPDYHKPLPPCR SVCERAKAGCSPLMRQYGFAPWPERMSCDRLPVLRDAEVLCDMDYNRSEATT APPRPFPAPKPTLPGPPGAPASGGGECPPAGGPFVCKCREPFVPIKESHPLYNKV RTGQVPNCAVPCYQPSFSADERTFATFWIGLWSVLCFISTSTTVATFLIDMERF RYPERPPIIFLSACYLCVSLGFLVRLVVGHASVACSREHNHIIHYETTGPALCTIVFL LVYFFGMASSIWWVILSLTWFLAAGMKWGNEAIAQYFHLAAWLIPSVKSIT ALALSSVDGDPVAGICYVGNQNLNSLRGFLVGLPLVLYLLVGTFLLAGFVSLFRI RSVIKQGGTKTDKLEKLMIRIGIFTLLYTPASIVVACYLEQHYRESWEAALTCA CPGHDTGQPRAKPEYWVLMKYFMCLVVGITSGVWWSGKTVESWRRFTSRC CCRPRRGHKSGGAMAAGDYPEASAALTGRTGPPGPAATYHKQVSLSHV
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	This antibody is reactive against transfected lysate in WB and as a detection antibody in ELISA.



Images

Western Blot: Frizzled-5 Antibody [H00007855-D01P] - Analysis of FZD5 expression in transfected 293T cell line by FZD5 polyclonal antibody. Lane 1: FZD5 transfected lysate(64.50 KDa). Lane 2: Non-transfected lysate.



Publications

Peterson YK, Nasarre P, Bonilla IV et al. Frizzled-5: a high affinity receptor for secreted frizzled-related protein-2 activation of nuclear factor of activated T-cells c3 signaling to promote angiogenesis. *Angiogenesis* 2017-08-24 [PMID: 28840375]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00007855-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

