

# Product Datasheet

## CDK2AP1 Antibody - Azide and BSA Free H00008099-B01P-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00008099-B01P-50ug](http://www.novusbio.com/H00008099-B01P-50ug)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00008099-B01P-50ug](http://www.novusbio.com/reviews/destination/H00008099-B01P-50ug)



**H00008099-B01P-50ug**

CDK2AP1 Antibody - Azide and BSA Free

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse CDK2AP1 Antibody - Azide and BSA Free (H00008099-B01P-50ug) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	8099
Gene Symbol	CDK2AP1
Species	Human
Immunogen	CDK2AP1 (AAH34717, 1 a.a. - 115 a.a.) full-length human protein. MSYKPNLAAHMPAAALNAAGSVHSPSTSMATSSQYRQLLSGYTQG TGNSQVPQSKYAELLAIEELGKEIRPTYAGSKSAMERLKRGIHARGLVRECLA ETERNARS
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/ Immunofluorescence




## Images

Western Blot: CDK2AP1 Antibody [H00008099-B01P-50ug] - Analysis of CDK2AP1 expression in transfected 293T cell line H00008099-T01) by CDK2AP1 MaxPab polyclonal antibody. Lane 1: CDK2AP1 transfected lysate(12.76 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: CDK2AP1 Antibody [H00008099-B01P-50ug] - Purified MaxPab antibody to CDK2AP1 on HepG2 cell. [antibody concentration 20 ug/ml]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00008099-B01P-50ug**

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBP1-40398	Recombinant Human CDK2AP1 His Protein

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00008099-B01P-50ug](http://www.novusbio.com/reviews/submit/H00008099-B01P-50ug)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

