

# Product Datasheet

## TAF1C Antibody (3E6) - Azide and BSA Free H00009013-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00009013-M02](http://www.novusbio.com/H00009013-M02)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00009013-M02](http://www.novusbio.com/reviews/destination/H00009013-M02)



**H00009013-M02**


TAF1C Antibody (3E6) - Azide and BSA Free

| Product Information         |   |
|-----------------------------|---|
| Unit Size                   | 0.1 mg  |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.  |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.  |
| Clonality                   | Monoclonal  |
| Clone                       | 3E6   |
| Preservative                | No Preservative   |
| Isotype                     | IgG2a Kappa   |
| Purity                      | IgG purified  |
| Buffer                      | In 1x PBS, pH 7.4   |
| Product Description         |   |
| Description                 | Quality control test: Antibody Reactive Against Recombinant Protein.  |
| Host                        | Mouse   |
| Gene ID                     | 9013  |
| Gene Symbol                 | TAF1C   |
| Species                     | Human   |
| Specificity/Sensitivity     | TAF1C - TATA box binding protein (TBP)-associated factor, RNA polymerase I, C, 110kDa (3E6)   |
| Immunogen                   | TAF1C (NP_005670.2, 761 a.a. ~ 869 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>SLSGHVDPS EDTSSPHSPEWPPADALPLPPTTPPSQELTPDACAQGV PSEQR<br>QMLRDYMAKLPPQRDTPGCATTPPHSQASSVRATRSQQHTPVLSSSQPLRKK<br>PRMGF |
| Notes                       | This product is produced by and distributed for Abnova, a company based in Taiwan.  |
| Product Application Details |   |
| Applications                | Western Blot, ELISA, Sandwich ELISA   |
| Recommended Dilutions       | Western Blot, ELISA, Sandwich ELISA   |
| Application Notes           | Antibody reactivity against transfected lysate and recombinant protein for WB. It has been used for ELISA.  |




## Images

Western Blot: TAF1C Antibody (3E6) [H00009013-M02] - Analysis of TAF1C expression in transfected 293T cell line by TAF1C monoclonal antibody (M02), clone 3E6. Lane 1: TAF1C transfected lysate(85.1 KDa). Lane 2: Non-transfected lysate.



Sandwich ELISA: TAF1C Antibody (3E6) [H00009013-M02] - Detection limit for recombinant GST tagged TAF1C is 1 ng/ml as a capture antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00009013-M02**

---

|                  |   |
|------------------|---|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK)                |
| NBP2-58841PEP    | TAF1C Recombinant Protein Antigen                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00009013-M02](http://www.novusbio.com/reviews/submit/H00009013-M02)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

