

# Product Datasheet

## MTMR4 293T Cell Transient Overexpression Lysate H00009110-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00009110-T01](http://www.novusbio.com/H00009110-T01)

Updated 8/21/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00009110-T01](http://www.novusbio.com/reviews/destination/H00009110-T01)



**H00009110-T01****MTMR4 293T Cell Transient Overexpression Lysate****Product Information**

|                      |                                                                                                           |
|----------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                                    |
| <b>Concentration</b> | mg/ml                                                                                                     |
| <b>Storage</b>       | Store at -80C. Avoid freeze-thaw cycles.                                                                  |
| <b>Buffer</b>        | 1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue) |

**Product Description**


|                                |                                                                                                                                                         |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Quality control test: Transient overexpression cell lysate was tested with Anti-MTMR4 antibody by Western Blots. <b>Plasmid:</b> pCMV-MTMR4 full-length |
| <b>Gene ID</b>                 | 9110                                                                                                                                                    |
| <b>Gene Symbol</b>             | MTMR4                                                                                                                                                   |
| <b>Species</b>                 | Human                                                                                                                                                   |
| <b>Specificity/Sensitivity</b> | MTMR4 293T Cell Transient Overexpression Lysate(Denatured)                                                                                              |
| <b>Notes</b>                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                      |
| <b>Lysate Type</b>             | 293T Cell Transient Overexpression                                                                                                                      |
| <b>Lysate Protein State</b>    | Denatured                                                                                                                                               |

**Product Application Details**


|                              |              |
|------------------------------|--------------|
| <b>Applications</b>          | Western Blot |
| <b>Recommended Dilutions</b> | Western Blot |

**Images**

Lane 1: MTMR4 transfected lysate (131.56 KDa) Lane 2: Non-transfected lysate.



MTMR4 transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00009110-T01](http://www.novusbio.com/reviews/submit/H00009110-T01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

