

# Product Datasheet

## SMC3 293T Cell Transient Overexpression Lysate H00009126-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00009126-T01](http://www.novusbio.com/H00009126-T01)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00009126-T01](http://www.novusbio.com/reviews/destination/H00009126-T01)



**H00009126-T01****SMC3 293T Cell Transient Overexpression Lysate****Product Information**

<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	mg/ml
<b>Storage</b>	Store at -80C. Avoid freeze-thaw cycles.
<b>Buffer</b>	1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue)

**Product Description**

<b>Description</b>	Quality control test: Transient overexpression cell lysate was tested with Anti-SMC3 antibody by Western Blots. <b>Plasmid:</b> pCMV-SMC3 full-length
<b>Gene ID</b>	9126
<b>Gene Symbol</b>	SMC3
<b>Species</b>	Human
<b>Specificity/Sensitivity</b>	SMC3 293T Cell Transient Overexpression Lysate(Denatured)
<b>Notes</b>	This product is produced by and distributed for Abnova, a company based in Taiwan.
<b>Lysate Type</b>	293T Cell Transient Overexpression
<b>Lysate Protein State</b>	Denatured

**Product Application Details**


<b>Applications</b>	Western Blot
<b>Recommended Dilutions</b>	Western Blot

**Images**

CSPG6 transfected lysate.



Lane 1: CSPG6 transfected lysate (141.5 KDa) Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00009126-T01**

---

H00009126-Q01-10ug	Recombinant Human SMC3 GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NB100-207	SMC3 Antibody
NB300-232	SCP3/SYCP3 Antibody - BSA Free

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00009126-T01](http://www.novusbio.com/reviews/submit/H00009126-T01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

