

Product Datasheet

PDLIM7 Antibody H00009260-B01

Unit Size: 0.05 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009260-B01

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009260-B01



H00009260-B01**PDLIM7 Antibody**

Product Information	
Unit Size	0.05 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	No additive


Product Description	
Host	Mouse
Gene ID	9260
Gene Symbol	PDLIM7
Species	Human
Specificity/Sensitivity	PDLIM7 - PDZ and LIM domain 7 (enigma),
Immunogen	PDLIM7 (NP_005442.2, 1 a.a. - 457 a.a.) full-length human protein. MDSFKVVL ^E GPAPWGFRLQGGKDFNVPLSISRLTPGGKAAQAGVAVGDWVLS IDGENAGSLTHIEAQNKIRACGERLSLGLSRAQPVQSKPQKASAPAADPPRYTF APSVSLNKTARPFGAPPPADSAPQQNGQPLRPLVPDASKQRLMENTEDWRPR PGTGQSR ^S FRILAHLTGTEFMQDPDEEHLKKSSQVPRTEAPAPASSTPQEPW PGPTAPSP ^T SRPPWAVDPFAFAERYAPDKTSTVLTRHSQPATPTPLQSRTSIVQ AAAGGVPGGGSNNGKTPVCHQCHKVIRGRYLVALGHAYHPEEFVCSQCGKVL EEGGFFEEKGAIFCPCYDVRYPSCAKCKKITGEIMHALKMTWHVHCFTCA ACKTPIRNRAFYMEEGVPYCERDYEKMF ^G TKCHGCD ^F KIDAGDRFLEALGFS WHDT ^C FVCAICQINLEGKTFY ^S KKDRPLCKSHAFSHV
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/ Immunofluorescence
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. It has also been used for immunofluorescence.




Images


Western Blot: PDLIM7 Antibody [H00009260-B01] - Analysis of PDLIM7 expression in transfected 293T cell line by PDLIM7 polyclonal antibody. Lane1:PDLIM7 transfected lysate(50.27 KDa). Lane2:Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: PDLIM7 Antibody [H00009260-B01] - Analysis of purified antibody to PDLIM7 on HeLa cell. Antibody concentration 10 ug/ml



Western Blot: PDLIM7 Antibody [H00009260-B01] - Western Blot analysis of PDLIM7 expression in transfected 293T cell line by PDLIM7 MaxPab polyclonal antibody (B01). Lane1: PDLIM7 transfected lysate (49.80 KDa). Lane2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009260-B01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

