

Product Datasheet

MD-1/LY86 Antibody H00009450-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009450-B01P

Updated 6/13/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009450-B01P



H00009450-B01P

MD-1/LY86 Antibody

Product Information

Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	IgG purified
Buffer	PBS (pH 7.4)

Product Description

Host	Mouse
Gene ID	9450
Gene Symbol	LY86
Species	Human
Specificity/Sensitivity	LY86 - lymphocyte antigen 86,
Immunogen	LY86 (NP_004262.1, 1 a.a. - 162 a.a.) full-length human protein. MKGFTATLFLWTLIFPSCSGGGGGKAWPTHVVCSDSGLEVLYQSCDPLQDFG FSVEKCSKQLKSNINIRFGIILREDIKELFLDLALMSQGSSVLNFSYPICEAALPKF SFCGRRKGGEQIYYAGPVNNPEFTIPQGEYQVLLLELYTEKRSTVACANATIMCS
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.


Product Application Details

Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500, Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	Antibody reactive against transfected lysate in western blot. May also be used for ELISA and immunohistochemistry (paraffin).




Images

Western Blot: MD-1/LY86 Antibody [H00009450-B01P] - Analysis of LY86 expression in transfected 293T cell line by LY86 polyclonal antibody. Lane 1: LY86 transfected lysate(17.82 KDa). Lane 2: Non-transfected lysate.



Immunohistochemistry-Paraffin: MD-1/LY86 Antibody [H00009450-B01P] - Analysis of purified antibody to LY86 on formalin-fixed paraffin-embedded human tonsil. (antibody concentration 3 ug/ml)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009450-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

