

Product Datasheet

BAG3 Antibody (5A8) **H00009531-M01A**

Unit Size: 0.2 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009531-M01A

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009531-M01A



H00009531-M01A


BAG3 Antibody (5A8)

Product Information	
Unit Size	0.2 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5A8
Preservative	No Preservative
Isotype	IgG2b Kappa
Purity	Unpurified
Buffer	Ascites
Product Description	
Host	Mouse
Gene ID	9531
Gene Symbol	BAG3
Species	Human
Immunogen	BAG3 (AAH14656.1, 1 a.a. ~ 575 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MSAATHSPMMQVASGNGDRDPLPPGWEIKIDPQTGWPFVVDHNSRTTTWND PRVPSEGPKETPSSANGPSREGSRLPPAREGHPVYPQLRPGYIPIVLHEGAE NRQVHPFHVYPQPGMQRFRTAAAAAPQRSQSPLRGMPETTQPDKQCGQVA AAAAAQPPASHGPERSQSPAASDCSSSSSSASLPSSGRSSLGSHQLPRGYISI PVIHEQNVTRPAAQPSFHQAQKTHYPAQQGEYQTHQPVYHKKIQGDDWEPRPL RAASPFRSSVQGASSREGSPARSSTPLHSPSPIRVHTVVDRPQQPMTTHRETA PVSQPENKPESKPGPVGPELPPGHIPIQVIRKEVDSKPVSKPPPPSEKVEVKV PPAPVPCPPSPGSPSAVPSSPKSVATEERAAPSTAPAEATPPKPGAEAPPKH PGVLKVEAILEKVQGLEQAVDNFEGKKTDDKKYLMIEEYLTKEALLDSVDPEGR ADVRQARRDGVKRVQTILEKLEQKAIDVPGQVQVYELQPSNLEADQPLQAIME MGAVAADKGKKNAGNAEDPHTETQQPEATAAATSNPSSMTDTPGNPAAP
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:100 - 1:2000, ELISA 1:100 - 1:2000



Images

Western Blot: BAG3 Antibody (5A8) [H00009531-M01A] - Analysis of BAG3 expression in transfected 293T cell line by BAG3 monoclonal antibody (M01A), clone 5A8. Lane 1: BAG3 transfected lysate (Predicted MW: 61.6 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00009531-M01A

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP1-72276-50ug	Recombinant Human BAG3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009531-M01A

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

