

Product Datasheet

Eukaryotic translation initiation factor 5B Antibody - Azide and BSA Free H00009669-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009669-B01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009669-B01P



H00009669-B01P

Eukaryotic translation initiation factor 5B Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description	
Host	Mouse
Gene ID	9669
Gene Symbol	EIF5B
Species	Human
Specificity/Sensitivity	EIF5B - eukaryotic translation initiation factor 5B,
Immunogen	EIF5B (AAH32639.1, 1 a.a. - 1220 a.a.) full-length human protein. MGKKQKNKSEDSTKDDIDLDAEIEGAGAAKEQEPQKSKGKKKKEKKKQDF DEDDILKELEELSLEAQGIKADRETVAVKPTENNEEEFTSKDKKKKGQKGGKQS FDDNDSEEDKDKSKSKTAKPKVEMYSGSDDDDDFNKLPKAKGKAQKSNK KWDGSEEDEDNSKKIKERSRINSSGESGDESDEFQSRKGQKKNQKNKPGPN IESGNEDDDASFIIKTVAQKKAERERKKRDEEKAKLRKLKEKEELETGKKD QSKQKESQRKFEEETVKS SVTVDTGVIPASEEKAETPTAAEDDNEGDKKKKDK KKKKGEKEEKEKEKKKGPSKATVKAMQEALAKLKEEEERQKREEERIKRLEE LEAKRKEEERLEQEKREKKQKEKERKERLKEGKLLTKSQREARARAEATLK LLQAQGVVPSKDSLPPKRPIYEDKKRKKIPQQLESKEVSESMELCAAVEVME QGVPEKEETPPPVEPEEEEDTEDAGLDDWEAMASDEETEKVEGNTVHIEVKE NPEEEEEEEEEEEEEEEEEEEEEEGESEGSEGDEEDEKVSDEKDSGKTLDK KPSKEMSSDSEYDSDDDRTEERAYDKAKRRIEKRRLEHSKNVNTEKLRAPIIC VLGHVDTGKTKILDKLRHHTHVQDGEAGGITQQIGATNVPLEAINEQTKMIKNFDR ENVRIIPGMLIIDTPGHESFSNLRNRGSSLCDIAILVVDIMHGLEPQTIESINLLKSK KCPFIVALKIDRLYDWKKSPDSVAATLKKQKKNTKDEFEERAKAIIVEFAQQG LNAALFYENKDPRTFVSLVPTSAHTGDGMGSLIYLLVELTQTMLSKRLAHCEEL RAQVMEVKALPGMGTTIDVILINGRLKEGDTIIVPGVEGPVITQIRGLLLPPPMKE LRVKNQYEKHKEVEAAQGVKILGKDLEKTLAGLPLL VAYKEDEIPVLKDELIHEL KQTLNAIKLEEKGVYVQASTLGSLALLEFLKTSEVPYAGINIGPVHKKDVMKAS VMLEHDPQYAVNLAFDVRIERDAQEMADSLGVRIFS AEIYHLFDAFTKYRQDY KKQKQEEFKHIAVFPCKIKILPQYIFNSRDPVIMGVTVEAGQVKQGTMPMCVPSK NFVDIGIVTSIEINHKQVDVAKKGQEVCKIEPIPGESPKMFGRHFEATDILVSKIS RQSIDALKDWFRDEM QKSDWQLIVELKKVFEII
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/ Immunofluorescence




Application Notes


Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.

Images


Western Blot: Eukaryotic translation initiation factor 5B Antibody [H00009669-B01P] - Analysis of EIF5B expression in human placenta.




Immunocytochemistry/Immunofluorescence: Eukaryotic translation initiation factor 5B Antibody [H00009669-B01P] - Analysis of purified antibody to EIF5B on HeLa cell. (antibody concentration 10 ug/ml)



Western Blot: Eukaryotic translation initiation factor 5B Antibody [H00009669-B01P] - Analysis of EIF5B expression in Jurkat.



Western Blot: Eukaryotic translation initiation factor 5B Antibody [H00009669-B01P] - Analysis of EIF5B expression in transfected 293T cell line by EIF5B polyclonal antibody. Lane 1: EIF5B transfected lysate (134.2 KDa). Lane 2: Non-transfected lysate.



Publications

Ma HY, Lu CY, Tsao KC et al. Association of EV71 3C polymorphisms with clinical severity. J Microbiol Immunol Infect 2017-06-28 [PMID: 28711442]



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00009669-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBP1-85120PEP	Eukaryotic translation initiation factor 5B Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009669-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

