

Product Datasheet

Hephaestin Antibody (2D3)

H00009843-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00009843-M01

Updated 10/13/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00009843-M01




H00009843-M01**Hephaestin Antibody (2D3)**


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2D3
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	9843
Gene Symbol	HEPH
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	HEPH (2D3)
Immunogen	HEPH (NP_620074, 315 a.a. ~ 425 a.a) partial recombinant protein with GST tag. TRGHHTDVANIFPATFVTAEMVPWEPGTWLISCQVNSHFRDGMQALYKVK SCSMAPPVDLLTGKVRQYFIEAHEIQWDYGPMGHGSGTKNLREPGSISDKFF QKSSSRI
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody reactive against cell lysate and recombinant protein for western blot. It has also been used for ELISA.

Images

Western Blot: Hephaestin Antibody (2D3) [H00009843-M01] - PH monoclonal antibody (M01), clone 2D3 Western Blot analysis of HEPH expression in K-562 (Cat # L009V1).



Sandwich ELISA: Hephaestin Antibody (2D3) [H00009843-M01] - Detection limit for recombinant GST tagged HEPH is approximately 0.3ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00009843-M01

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981	Mouse IgG2a Isotype Control (M2AK)
H00009843-Q01-10ug	Recombinant Human Hephaestin Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00009843-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

