

Product Datasheet

Optineurin Antibody - Azide and BSA Free H00010133-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010133-B01P

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010133-B01P



H00010133-B01P

Optineurin Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Mouse Optineurin Antibody - Azide and BSA Free (H00010133-B01P) is a polyclonal antibody validated for use in WB, ICC/IF and Simple Western. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	10133
Gene Symbol	OPTN
Species	Human
Specificity/Sensitivity	OPTN - optineurin,
Immunogen	OPTN (NP_001008212.1, 1 a.a. - 577 a.a.) full-length human protein. MSHQPLSCLTEKEDSPSESTGNGPPHLAHPNLDFTPEELLQQMKELLTENHQ LKEAMKLNNQAMKGRFEELSAWTEKQKEERQFFEIQSKEAKERLMALSHENE KLKEELGKLKGKSERSSDPTDDSRLPRAEAEQEKDQLRTQVVRLQAEKADLL GIVSELQLKLNSSGSSSEDSFVEIRMAEGEAEGSVKEIKHSPGPTRTVSTGTALS KYRSRSADGAKNYFEHEELTVSQQLLCLREGNQKVERLEVALKEAKERVSDFE KKTSNRSEIETQTEGSTEKENDEEKGPETVGSEVEALNLQVTSLFKELQEAHTK LSEAELMKKRLQEKCQALERKNSAIPSELNEKQELVYTNKKLELQVESMLSEIK MEQAKTEDEKSKLTVLQMTHNKLQEHNNALKTIEELTRKESEKVDRAVLKELS EKLELAEKALASKQLQMDKQTIKQEEDELETMTILRAQMEVYCSDFAERA AREKIHEEKEQLALQLAVLLKENDAFEDGGRQSLMEMQSRHGARTSDSDQQA YLVQRGAEDRDWRQQRNIPIHSCPCKGCVLPDIDTLQIHVMDCII
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Simple Western, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Simple Western, Immunocytochemistry/ Immunofluorescence 1:10 - 1:500
Application Notes	This Optineurin Antibody is validated for Simple Western from a verified customer review. See Simple Western Antibody Database for Simple Western validation: Tested in THP-1 macrophage lysate 0.2 mg/mL, 0.6 mg/mL and 1.2 mg/mL, separated by Size, antibody dilution of 1:25




Images

Western Blot: Optineurin Antibody [H00010133-B01P] - Analysis of OPTN expression in human placenta.




Immunocytochemistry/Immunofluorescence: Optineurin Antibody [H00010133-B01P] - Analysis of purified antibody to OPTN on HeLa cell. Antibody at 10 ug/mL.



Western Blot: Optineurin Antibody [H00010133-B01P] - Analysis of OPTN expression in transfected 293T cell line by OPTN polyclonal antibody. Lane 1: OPTN transfected lysate (63.47 kDa). Lane 2: Non-transfected lysate.



Simple Western: Optineurin Antibody [H00010133-B01P] - This antibody (1:25 dilution) was used to probe THP-1 macrophage lysate (concentrations shown) by Simple Western. Simple Western image submitted by a verified customer review.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00010133-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBP1-84682PEP	Optineurin Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010133-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

