

Product Datasheet

A33/GPA33 Antibody H00010223-B01

Unit Size: 0.05 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010223-B01

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010223-B01



H00010223-B01

A33/GPA33 Antibody

Product Information	
Unit Size	0.05 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Unpurified
Buffer	No additive


Product Description	
Host	Mouse
Gene ID	10223
Gene Symbol	GPA33
Species	Human
Specificity/Sensitivity	GPA33 - glycoprotein A33 (transmembrane),
Immunogen	GPA33 (NP_005805.1, 1 a.a. - 319 a.a.) full-length human protein. MVGKMWPVLWTLCAVRVTVDASVETPQDVLRSQGKSVTLPTCTYHTSTSSR EGLIQWDKLLLTHTERVVIWPFNSKNYIHGELYKNRVSISSNNAEQSDASITIDQLT MADNGTYECSVSLMSDLEGNTKSRVLLVLPVPSKPECGIEGETIIGNNIQLTC QSKEGSPTPQYSWKRYNILNQEQLAQPASGQPVSLKNISTDTSGYYICTSSN EEGTQFCNITVAVRSPSMNVALYVGVVVAALIIIGIIIYCCCRGKDDNTEDK EDARPNREAYEEPPEQLRELSREREEDDYRQEEQRSTGRES PDHLDQ
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500, Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	Antibody reactive against transfected lysate and tissue lysate for Western Blot. Has also been used for immunohistochemistry (paraffin) and ELISA.




Images


Western Blot: A33/GPA33 Antibody [H00010223-B01] - Analysis of GPA33 expression in human colon.



Immunohistochemistry-Paraffin: A33/GPA33 Antibody [H00010223-B01] - Analysis of antibody to GPA33 on formalin-fixed paraffin-embedded human small Intestine. Antibody concentration 3 ug/ml



Western Blot: A33/GPA33 Antibody [H00010223-B01] - Analysis of GPA33 expression in transfected 293T cell line. Lane 1: GPA33 transfected lysate (35.09 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010223-B01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

