

Product Datasheet

UBAC1 Antibody H00010422-B01P

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010422-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010422-B01P



H00010422-B01P

UBAC1 Antibody

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	10422
Gene Symbol	UBAC1
Species	Human
Reactivity Notes	This product is reactive against Human.
Immunogen	UBAC1 (NP_057256.2, 1 a.a. - 405 a.a.) full-length human protein. MFVQEEKIFAGKVLRLHICASDGAEWLEEATEDTSVEKLERCLKHCAHGSLE DPKSITHHKKLIHAASERVLSDARTILEENIQDQDVLLLIKKRAPSPKPMADVSAE EKKKQDQKAPDKEAILRATANLPSYNMDRAAVQTNMRDFQTELRKILVSLIEVA QKLLALNPDAVELFKKANAMLDEDEDERVDAAALRQLTEMGFPENRATKALQL NHMSVPQAMEWLIEHAEDPTIDTPLPGQAPPEAEGATAAASEAAAGASATDEE ARDELTEIFKIRRKREFRADARAVISLMEMGFDEKEVIDALRVNNNQNAACE WLLGDRKPSPEELDKGIDPDSPLFQAILDNPVVQLGLTNPKTLLAFEDMLENPL NSTQWMNDPETGPVMLQISRIFQTLNRT
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/Immunofluorescence
Application Notes	This antibody is useful for Western Blot




Images


Western Blot: UBAC1 Antibody [H00010422-B01P] - Analysis of UBADC1 expression in human liver.



Immunocytochemistry/Immunofluorescence: UBAC1 Antibody [H00010422-B01P] - Analysis of purified antibody to UBAC1 on HeLa cell. (antibody concentration 10 ug/ml)



Western Blot: UBAC1 Antibody [H00010422-B01P] - Analysis of UBAC1 expression in transfected 293T cell line by UBAC1 polyclonal antibody. Lane 1: UBADC1 transfected lysate(44.55 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010422-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

