

Product Datasheet

hnRNP-Q Antibody (2E10)

H00010492-M05

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00010492-M05

Updated 7/6/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00010492-M05



H00010492-M05


hnRNP-Q Antibody (2E10)

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2E10
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	Protein A purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	10492
Gene Symbol	SYNCRIP
Species	Human
Specificity/Sensitivity	Reacts with synaptotagmin binding, cytoplasmic RNA interacting protein.
Immunogen	SYNCRIP (AAH32643, 1 a.a. ~ 527 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MATEHVNGNGTEEPMDTTSAVIHSENFQTLDDAGLPQKVAEKLDEIYVAGLVA HSDLDERAIEALKEFNEDGALAVLQQFKDSLHVNKSAFLCGVMKTYRQRE KQGTKVADSSKGPDEAKIKALLERTGYTLDDVTTGQRKYGGPPPSVYSGQQP SVGTEIFVKGIPRDLFEDELVPLFEKAGPIWDLRLMMDPLTGLNRYAFVTFCT KEAAQEAVKLYNNHEIRSGKHIGVCISVANNRLFVGSIPKSKTKEQILEEFSKVT EGLTDVILYHQPDDKKKNRGFCFLEYEDHKTAQVKVLFVRNLANTVTEEILEK AFSQFGKLERVKKLKDYAFIHFDERDGTVKAMEEMNGKDLEGENIEIVFAKPPD QKRKERKAQRQAANKMYDDYYYYGPPHMPPTTRGRGRGGRGGYGYPPDY YGYEDYYDYGYDYHNYRGGYEDPYGYEDFQVGARGRGGRGARGAAPS GRGAAPPRGRAGYSQRGGPGSARGVARGGAQQQRGRGQKGV EAGPD LLQ
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	ELISA, Sandwich ELISA
Recommended Dilutions	ELISA, Sandwich ELISA
Application Notes	This antibody is reactive against recombinant protein in ELISA.



Images

Sandwich ELISA: hnRNP-Q Antibody (2E10) [H00010492-M05] -
Detection limit for recombinant GST tagged SYNCRIP is 0.1 ng/ml as a
capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00010492-M05

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
H00010492-P02-2ug	Recombinant Human hnRNP-Q GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00010492-M05

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

