

Product Datasheet

Septin-8 Antibody (1G3-7D1)

H00023176-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00023176-M01

Updated 5/13/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00023176-M01




H00023176-M01**Septin-8 Antibody (1G3-7D1)**

Product Information	
Unit Size	0.1 mg
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1G3-7D1
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Ascites
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	23176
Gene Symbol	SEPTIN8
Species	Human
Immunogen	SEPT8 (AAH01329, 1 a.a. ~ 258 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MKKLDISKVNIPIIAKADTISKSELHKFKIKIMGELVSNGVQIYQFPTDDEAVAEINA VMNAHLPPFAVVGSTEEVKVGNKLVRARQYPWGVVQVENENHCDFVKLREMLI RVNMEDLREQTHSRHYELYYRCKLEEMGFQSDGDSQPFSLQETYEAKRKEF LSELQRKEEEMRQMFVNKVKETELKEKERELHEKFEHLKRVHQEEKRKVEE KRRELEEETNAFNRRKAAVEALQSQUALHATSQQPLRKDKDKKN
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Sandwich ELISA
Recommended Dilutions	Western Blot 1:100 - 1:2000, ELISA 1:100 - 1:2000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:10 - 1:500, Sandwich ELISA




Images


Immunocytochemistry/Immunofluorescence: Septin-8 Antibody (1G3-7D1) [H00023176-M01] - Analysis of monoclonal antibody to SEPT8 on HeLa cell . Antibody concentration 10 ug/ml.



Immunohistochemistry-Paraffin: Septin-8 Antibody (1G3-7D1) [H00023176-M01] - Analysis of monoclonal antibody to SEPT8 on formalin-fixed paraffin-embedded human small Intestine. Antibody concentration 3 ug/ml.



Sandwich ELISA: Septin-8 Antibody (1G3-7D1) [H00023176-M01] - Detection limit for recombinant GST tagged SEPT8 is 3 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00023176-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

