

Product Datasheet

MKRN2 Antibody **H00023609-B01P**

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00023609-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00023609-B01P



H00023609-B01P**MKRN2 Antibody**

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	23609
Gene Symbol	MKRN2
Species	Human
Specificity/Sensitivity	MKRN2 - makorin, ring finger protein, 2,
Immunogen	MKRN2 (NP_054879.3, 1 a.a. - 416 a.a.) full-length human protein. MSTKQITCRYFMHGVCREGSQCLFSDLANSPSTICKYYQKGYCAYGTRCR YDHTRPSSAAAGGAVGTMAHSVPSPAFHSPHPPSEVTASIVKTNSHEPGKREK RTLVLDRNLSGMAERKTQPSMVSNPGSCSDPQPSPKPHSYLDAIRSGLD DVEASSSYSNEQQLCPYAAAGECRFGDACVYLHGEVCEICRLQVLHPFDPEQ RKAHEKICMLTFEHEMEKAFQASQDKVCSICMEVILEKASASERRFGILSNC NHTYCLSCIRQWRCQKQFENPIIKSCPECRVISEFVIPSVMYVWEDQNKKNELIEA FKQGMGKKACKYFEQKGKTCFPGSKCLYRHAYPDGRLAEPEKPRKQLSSQG TVRFFNSVRLWDFIENRESRHVPNNEDVDMTELGDLMHLSGVESSEP
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500
Application Notes	Antibody reactive against tissue lysate and transfected lysate in western blot. May also be used for ELISA.




Images

Western Blot: MKRN2 Antibody [H00023609-B01P] - Analysis of MKRN2 expression in human liver.



Western Blot: MKRN2 Antibody [H00023609-B01P] - Analysis of MKRN2 expression in transfected 293T cell line by MKRN2 polyclonal antibody. Lane 1: MKRN2 transfected lysate(45.76 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00023609-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

