

Product Datasheet

PRPF31 Antibody H00026121-B02P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00026121-B02P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00026121-B02P



H00026121-B02P

PRPF31 Antibody


Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)

Product Description	
Host	Mouse
Gene ID	26121
Gene Symbol	PRPF31
Species	Human
Specificity/Sensitivity	This product is specific for Human PRPF31 purified MaxPab mouse polyclonal antibody (B02P) [Gene ID: 26121].
Immunogen	PRPF31 (NP_056444.2, 1 a.a. - 499 a.a.) full-length human protein. MSLADELLADLEEAEEEEEGGSYGEEEEPAIEDVQEETQLDLSGDSVKTIACL WDSKMF AEIMMKIEEYISKQAKASEVMGPVEAAPEYRVIVDANNLTVEIENELNII HKFIRDKYSKRFP ELES LVPNALDYIRTVKELGNSLDKCKNNENLQQILT NATIM VVSVTASTTQGGQLSEEELEERLEEACDMALELNASKHRIY EYVESRMSFIAPNL SIIIGASTAAKIMGVAGGLTNLSKMPACNIMLLGAQRKTLSGFSSTSVLPHTGYIY HSDIVQSLPPDLRRKAARLVAAKCTLAARVDSFH ESTEGKVG YELKDEIERKFD KWQEP PPVKQVKPLPAPLDGQRKRGRRYRKMKERLGLTEIRKQANRMSF GEIEEDAYQEDLG FSLGHLGKSGSGRVRQTQVNEATKARISKTLQRTLQKQSV VYGGKSTIRDRSSGTASSVAFTPLQGLEIVNPQAAEKKVAEANQKYFSSMAEFL KVKGEKSGLMST
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	Mouse polyclonal antibody raised against a full-length human PRPF31 protein.

Images

Western Blot: PRPF31 Antibody [H00026121-B02P] - Analysis of PRPF31 expression in transfected 293T cell line by PRPF31 polyclonal antibody. Lane 1: PRPF31 transfected lysate(55.5 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00026121-B02P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

