

# Product Datasheet

## TAF5L Antibody - Azide and BSA Free H00027097-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00027097-D01P](http://www.novusbio.com/H00027097-D01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00027097-D01P](http://www.novusbio.com/reviews/destination/H00027097-D01P)



**H00027097-D01P**


TAF5L Antibody - Azide and BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 mg                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                     |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                                         |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Preservative                | No Preservative                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Purity                      | Protein A purified                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Buffer                      | PBS (pH 7.4)                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Description                 | Novus Biologicals Rabbit TAF5L Antibody - Azide and BSA Free (H00027097-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                                                                  |
| Host                        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Gene ID                     | 27097                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Gene Symbol                 | TAF5L                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Specificity/Sensitivity     | TAF5L - TAF5-like RNA polymerase II, p300/CBP-associated factor (PCAF)-associated factor, 65kDa,                                                                                                                                                                                                                                                                                                                                     |
| Immunogen                   | TAF5L (NP_001020418.1, 1 a.a. - 325 a.a.) full-length human protein.<br>MKRVRTEQIQMAVSCYLKRRQYVDSGPKQGLRLSQTAEEEMAANLTVQSES<br>GCANIVSAAPCQAEPQQYEVQFGRLRNFLTSDSQHSHEVMPLLYPLFVYLHL<br>NLVQNSPKSTVESFYSRFHGMFLQNASQKDVIEQLQTTQTIQDILSNFKLRAFL<br>DNKYVVRLQEDSYNYLIRYLQSDNNTALCKVLTLLHIHLDVQPAKRTDYQLYASG<br>SSSRSENNGLEPPDMPSPILQNEAALEVLQESIKRVKDGPPSLTTICFYAFYNTE<br>QLLNTAEISPDSKLLAAGFDNSCIKLWSLRSKKLKSEPHQVDVSRIHLACDILEE<br>EV |
| Notes                       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                                   |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Recommended Dilutions       | Western Blot 1:500                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Application Notes           | Antibody reactive against recombinant protein with GST tag on ELISA and western blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.                                                                                                                                                                                                                                                    |




## Images

Western Blot: TAF5L Antibody [H00027097-D01P] - Analysis of TAF5L expression in HeLa.



Western Blot: TAF5L Antibody [H00027097-D01P] - Analysis of TAF5L expression in transfected 293T cell line by TAF5L polyclonal antibody. Lane 1: TAF5L transfected lysate(37.00 KDa). Lane 2: Non-transfected lysate.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to H00027097-D01P**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00027097-P01-10ug | Recombinant Human TAF5L GST (N-Term) Protein        |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00027097-D01P](http://www.novusbio.com/reviews/submit/H00027097-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

