

Product Datasheet

PR48 Antibody H00028227-B02P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00028227-B02P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00028227-B02P



H00028227-B02P

PR48 Antibody

Product Information

Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description

Host	Mouse
Gene ID	28227
Gene Symbol	PPP2R3B
Species	Human
Specificity/Sensitivity	PPP2R3B - protein phosphatase 2 (formerly 2A), regulatory subunit B', beta,
Immunogen	PPP2R3B (NP_955358.1, 1 a.a. - 176 a.a.) full-length human protein. MDLDGDGALSMFELEYFYEEQCRRRLDSMAIEALPFQDCLCQMLDLVKPRTEGK ITLQDLKRCKLANVFFDTFFNIEKYLDHEQKEQISLLRDGDSGGPELSDWEKYA AEEYDILVAEETAGEPWEDGF EAELSPVEQKLSALRSPLAQRPF FEAPSPLGAV DLYEYACGDEDLEPL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.


Product Application Details

Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/Immunofluorescence
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. It has also been used for immunofluorescence.




Images

Western Blot: PR48 Antibody [H00028227-B02P] - Analysis of PPP2R3B expression in transfected 293T cell line by PPP2R3B polyclonal antibody. Lane 1: PPP2R3B transfected lysate(19.36 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: PR48 Antibody [H00028227-B02P] - Analysis of purified antibody to PPP2R3B on HeLa cell. (antibody concentration 10 ug/ml)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00028227-B02P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

