

Product Datasheet

SART2 Antibody (6D4) H00029940-M03

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00029940-M03

Updated 10/13/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00029940-M03




H00029940-M03**SART2 Antibody (6D4)**


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6D4
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	29940
Gene Symbol	DSE
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	SART2 - squamous cell carcinoma antigen recognized by T cells 2
Immunogen	SART2 (NP_037484, 574 a.a. ~ 674 a.a) partial recombinant protein with GST tag.LGEESPLETEAASFFHNVDVPFEETVVDGVHGFIRQRDGLYKMYWMDDTG YSEKATFASVTYPRGYNGTNYVNVMTMLRSPITRAAYLFIGPSIDVQS
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:500, ELISA
Application Notes	Antibody reactive against recombinant protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.

Images

Western Blot: SART2 Antibody (6D4) [H00029940-M03] - Analysis of DSE expression in transfected 293T cell line by DSE monoclonal antibody (M03), clone 6D4. Lane 1: DSE transfected lysate (Predicted MW: 109.8 KDa). Lane 2: Non-transfected lysate.



ELISA: SART2 Antibody (6D4) [H00029940-M03] - Detection limit for recombinant GST tagged SART2 is approximately 0.3ng/ml as a capture antibody.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00029940-M03

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43318-0.5mg	Mouse IgG2a Kappa Light Chain Isotype Control (M2a)
H00029940-P01-2ug	Recombinant Human SART2 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00029940-M03

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

