

Product Datasheet

Neurogenin-3 Antibody (4B7)

H00050674-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00050674-M01

Updated 2/1/2019 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00050674-M01




H00050674-M01**Neurogenin-3 Antibody (4B7)**


Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4B7
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	50674
Gene Symbol	NEUROG3
Species	Human
Reactivity Notes	Human. Other species not tested.
Specificity/Sensitivity	NEUROG3 - neurogenin 3
Immunogen	NEUROG3 (NP_066279, 1 a.a. ~ 101 a.a) partial recombinant protein with GST tag.MTPQPSGAPTQVTRETERSFPRASEDEVTCPTSAPPSPTRTRGNCAEAE EGGCRGAPRKLRRARRGGRSRPKSELALSKQRRSRRKKANDRERNRMHNLN
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Application Notes	Antibody Reactive against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.

Images

Immunocytochemistry/Immunofluorescence: Neurogenin-3 Antibody (4B7) [H00050674-M01] - Analysis of monoclonal antibody to NEUROG3 on HeLa cell . Antibody concentration 10 ug/ml.



Immunoprecipitation: Neurogenin-3 Antibody (4B7) [H00050674-M01] - Analysis of NEUROG3 transfected lysate using anti-NEUROG3 monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with NEUROG3 MaxPab rabbit polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00050674-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

