

# Product Datasheet

## VE-Statin/EGFL7 Antibody

### H00051162-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00051162-B01P](http://www.novusbio.com/H00051162-B01P)

Updated 1/12/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00051162-B01P](http://www.novusbio.com/reviews/destination/H00051162-B01P)



**H00051162-B01P****VE-Statin/EGFL7 Antibody****Product Information**

<b>Unit Size</b>	0.05 mg
<b>Concentration</b>	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	No Preservative
<b>Isotype</b>	IgG
<b>Purity</b>	Protein G purified
<b>Buffer</b>	PBS (pH 7.4)

**Product Description**


<b>Host</b>	Mouse
<b>Gene ID</b>	51162
<b>Gene Symbol</b>	EGFL7
<b>Species</b>	Human
<b>Specificity/Sensitivity</b>	This product is specific for Human EGFL7 purified MaxPab mouse polyclonal antibody (B01P) [Gene ID: 51162].
<b>Immunogen</b>	EGFL7 (NP_057299.1, 1 a.a. - 273 a.a.) full-length human protein. MRGSQEVLLMWLLVLAVGGTEHAYRPGRRVCAVRAHGDVPVSESFVQRVYQP FLTTCDGHRACSTYRTIYRTAYRRSPGLAPARPRYACCPGWKRTSGLPGACG AAICQPPCRNGGSCVQPGRCRCPAGWRGDTQCSDVDECSARRGGCPQRCV NTAGSYWCQCWEGHSLADGTLCVPKGGPPRVAPNPTGVDSAMKEEVQRLQ SRVDLLEEKLQLVLAPLHSLASQALEHGLPDPGSLLVHSFQQLGRIDSLSEQISF LEEQLGSCSCKKDS
<b>Notes</b>	This product is produced by and distributed for Abnova, a company based in Taiwan.

**Product Application Details**

<b>Applications</b>	Western Blot
<b>Recommended Dilutions</b>	Western Blot
<b>Application Notes</b>	Mouse polyclonal antibody raised against a full-length human EGFL7 protein.

**Images**

Western Blot: VE-Statin/EGFL7 Antibody [H00051162-B01P] - Analysis of EGFL7 expression in transfected 293T cell line by EGFL7 polyclonal antibody. Lane 1: EGFL7 transfected lysate(30.03 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis.  
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00051162-B01P](http://www.novusbio.com/reviews/submit/H00051162-B01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

