

Product Datasheet

ETS1 associated protein II Antibody H00051567-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00051567-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00051567-B01P



H00051567-B01P

ETS1 associated protein II Antibody

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	51567
Gene Symbol	TDP2
Species	Human, Mouse
Specificity/Sensitivity	Reacts with TRAF and TNF receptor associated protein.
Immunogen	TTRAP (NP_057698.2, 1 a.a. - 362 a.a.) full-length human protein. MELGSCLEGGREAAEEEGEPEVKKRRLLCVEFASVASCDAAVAQCFLAENDW EMERALNSYFEPPEVESALERRPETISEPKTYVDLTNEETDSTTSKISPSEDQ QENGSFMFLITWNIDGLDLNLSERARGVCSYLALYSPDVIFLQEVIPPYYSYLK KRSSNYEITGHEEGYFTAIMLKKSRVKLKSQEIIPFPSTKMMRNLLCVHVNVS NELCLMTSHLESTRGHAAERMNQLKMLVKKMQEAPESATVIFAGDTNLRDREV TRCGGLPNNIVDVWEFLGKPKHCQYTWDTQMNSNLGITAACKLRFDRIFFRAA AEEGHIIPRSLDLLGLEKLDCGRFSDHWGLLCNLDIIL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/Immunofluorescence
Application Notes	This antibody is reactive against transfected lysate in western blot, immunofluorescence, and as a detection antibody in ELISA.




Images

Western Blot: ETS1 associated protein II Antibody [H00051567-B01P] - Analysis of TTRAP expression in transfected 293T cell line by TTRAP polyclonal antibody. Lane 1: TTRAP transfected lysate(39.82 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: ETS1 associated protein II Antibody [H00051567-B01P] - Analysis of purified antibody to TTRAP on Hs 181.Tes cell. (antibody concentration 10 ug/ml)



Publications

Thomson G, Watson A, Caldecott K et al. Generation of assays and antibodies to facilitate the study of human 5'-tyrosyl DNA phosphodiesterase. *Anal Biochem.* 2013-02-04 [PMID: 23416181]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00051567-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

