

Product Datasheet

CPSF73 Antibody H00051692-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00051692-D01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00051692-D01P



H00051692-D01P**CPSF73 Antibody**

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Rabbit
Gene ID	51692
Gene Symbol	CPSF3
Species	Human
Reactivity Notes	This product is reactive against Human.
Immunogen	CPSF3 (NP_057291.1, 1 a.a. - 684 a.a.) full-length human protein. MSAIPAEESDQLLIRPLGAGQEVGRSCIILEFKGRKIMLDCGIHPGLEGM DALPYI DLIDPAEIDLLLSHFHLDHCGALPWFLQKTSFKGRTFMTHATKAIYRWLLSDYV KVSNISADDMLYTETDLEESMDKIETINFHEVKEVAGIKFWCYHAGHV LGAAMF MIEIAGVKLLYTGDFSRQEDRHLMAAEIPNIKPDILIIESTYGTHIHEKREEREARF CNTVHDIIVNRGGRGLIPVFALGRAQELLILDEYWQNHPELHDIPYYASSLAKK CMAVYQTYVNAMNDKIRKQININNPVFVKHISNLKSMDHFDDIGPSVVMASPGM MQSGLSRELFEWCTDKRNGVIIAGYCVEGTLAKHIMSEPEEITMSGQKLPLK MSVDYISFSAHTDYQQTSEFIRALKPPHVILVHGEQNEMARLKAALIREYEDND EVHIEVHNPRNTEAVTLNFRGEKLAKVMGFLADKKPEQGQQRVSGILVKRNFNY HILSPCDLSNYTDLAMSTVKQTQAIPYTGPFNLLCYQLQKLTGDVEELEIQEKPA LKVFKNITVIQEPGMVVLEWLANPSNDMYADTVTTVILEVQSNPKIRKGAVQKV SKKLEMHVYSKRLEIMLQDIFGEDCVSVKDDSILSVTVDGKTANLNLETRTVEC EEGSEDDESLREMV ELAAQRLYEALTPVH
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot
Application Notes	This antibody is useful for Western Blot




Images

Western Blot: CPSF73 Antibody [H00051692-D01P] - Analysis of CPSF3 expression in HeLa.



Western Blot: CPSF73 Antibody [H00051692-D01P] - Analysis of CPSF3 expression in transfected 293T cell line by CPSF3 polyclonal antibody. Lane 1: CPSF3 transfected lysate (77.50 kDa). Lane 2: Non-transfected lysate.



Publications

Johansson C, Somberg M, Li X et al. HPV-16 E2 contributes to induction of HPV-16 late gene expression by inhibiting early polyadenylation. EMBO J. 2012-05-22 [PMID: 22617423]

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00051692-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

