

# Product Datasheet

## TULA/STS-2 Antibody H00053347-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00053347-D01P](http://www.novusbio.com/H00053347-D01P)

Updated 1/12/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00053347-D01P](http://www.novusbio.com/reviews/destination/H00053347-D01P)



**H00053347-D01P**

TULA/STS-2 Antibody

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 mg                                                                                                           |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | No Preservative                                                                                                  |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Protein A purified                                                                                               |
| <b>Buffer</b>        | PBS (pH 7.4)                                                                                                     |

**Product Description**

|                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>                    | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Gene ID</b>                 | 53347                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Gene Symbol</b>             | UBASH3A                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Species</b>                 | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Specificity/Sensitivity</b> | UBASH3A - ubiquitin associated and SH3 domain containing, A,                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen</b>               | UBASH3A (AAH28138.1, 1 a.a. - 451 a.a.) full-length human protein.<br>MAAGETQLYAKVSNKLKGRSSPSLLEPLLAMGFPVHTALKALAATGRKTAEEA<br>LAWLHDHCNDPSLDDPIPQEYALFLCPTGPLLEKLQEFWRESKRQCAKNRAHE<br>VFPHVTLCDDFTCEDQKVECLYEALKRAGDRLLGSFPTAVPLALHSSISYLGFF<br>VSGSPADVIREFAMTFATEASLLADCSVKPCTKQLHLTLAHKFYPHHQRTLEQL<br>ARAIPLGHSCQWTAALYSRDMRFVHYQTLRALFQYKPQNVDELTLSPGDYIFV<br>DPTQQDEASEGWVIGISQRTGCRGFLPENYTDRAESDWTWKHRMYTFSLAT<br>DLNSRKDGEASSRCSGEFLPQTARSLSSLQALQATVARKSVLVVRHGERVDQI<br>FGKAWLQQCSTPDGNYRDPDLNFPCLPRRSRGIKDFENDPPLSSCGIFQSRI<br>AGMFEDCLVGKVTITATLIMASRHPA |
| <b>Notes</b>                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |


**Product Application Details**

|                              |                                                                                                                                                                                   |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot                                                                                                                                                                      |
| <b>Recommended Dilutions</b> | Western Blot 1:500                                                                                                                                                                |
| <b>Application Notes</b>     | Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. |



## Images

Western Blot: TULA/STS-2 Antibody [H00053347-D01P] - Analysis of UBASH3A expression in transfected 293T cell line by UBASH3A polyclonal antibody. Lane 1: UBASH3A transfected lysate(50.40 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis.  
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00053347-D01P](http://www.novusbio.com/reviews/submit/H00053347-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

