

Product Datasheet

DNA Polymerase epsilon subunit 3 293T Cell Transient Overexpression Lysate H00054107-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00054107-T01

Updated 8/23/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00054107-T01




H00054107-T01**DNA Polymerase epsilon subunit 3 293T Cell Transient Overexpression Lysate**


Product Information	
Unit Size	0.1 ml
Concentration	mg/ml
Storage	Store at -80C. Avoid freeze-thaw cycles.
Buffer	1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue)
Product Description	
Description	Plasmid: pCMV-POLE3 full-length
Gene ID	54107
Gene Symbol	POLE3
Species	Human
Specificity/Sensitivity	DNA Polymerase epsilon subunit 3 293T Cell Transient Overexpression Lysate (Denatured)
Notes	<p>Quality control test: Transient overexpression cell lysate was tested with Anti-POLE3 antibody by Western Blots.</p> <p>This product is produced by and distributed for Abnova, a company based in Taiwan.</p>
Lysate Type	293T Cell Transient Overexpression
Lysate Protein State	Denatured
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot

Images

POLE3 transfected lysate



Lane 1: POLE3 transfected lysate (16.28 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00054107-T01

NBP2-23341	Recombinant Human DNA Polymerase epsilon subunit 3 His Protein
NB100-56708	Caspase-3 Antibody (31A1067) - (Pro and Active)
NBP3-03281-100ul	DNA Polymerase epsilon subunit 3 Antibody
DRT100	TNF RI/TNFRSF1A [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00054107-T01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

