

Product Datasheet

PTPIP51 293T Cell Transient Overexpression Lysate H00055177-T03

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00055177-T03

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00055177-T03



H00055177-T03**PTPIP51 293T Cell Transient Overexpression Lysate**

Product Information	
Unit Size	0.1 ml
Concentration	mg/ml
Storage	Store at -80C. Avoid freeze-thaw cycles.
Buffer	1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue)

Product Description	
Description	Quality control test: Transient overexpression cell lysate was tested with Anti-FAM82C antibody by Western Blots. Plasmid: pCMV-FAM82C full-length
Gene ID	55177
Gene Symbol	RMDN3
Species	Human
Specificity/Sensitivity	FAM82C 293T Cell Transient Overexpression Lysate(Denatured)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Lysate Type	293T Cell Transient Overexpression
Lysate Protein State	Denatured

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00055177-T03

H00055177-P01-10ug	Recombinant Human PTPN22 GST (N-Term) Protein
236-EG-200	EGF [Unconjugated]
NBP1-84738	PTN22 Antibody
AF231	EGFR Antibody [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00055177-T03

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

