

# Product Datasheet

## CSTP1 Antibody - Azide and BSA Free H00055313-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00055313-B01P](http://www.novusbio.com/H00055313-B01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00055313-B01P](http://www.novusbio.com/reviews/destination/H00055313-B01P)



**H00055313-B01P**


CSTP1 Antibody - Azide and BSA Free

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.05 mg                                                                                                                                                                                                                                                                                                                                                                                                            |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                   |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                                                                                                       |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                         |
| Preservative                | No Preservative                                                                                                                                                                                                                                                                                                                                                                                                    |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                                                                                                |
| Purity                      | Protein A purified                                                                                                                                                                                                                                                                                                                                                                                                 |
| Buffer                      | PBS (pH 7.4)                                                                                                                                                                                                                                                                                                                                                                                                       |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Description                 | Novus Biologicals Mouse CSTP1 Antibody - Azide and BSA Free (H00055313-B01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene ID                     | 55313                                                                                                                                                                                                                                                                                                                                                                                                              |
| Gene Symbol                 | CPPED1                                                                                                                                                                                                                                                                                                                                                                                                             |
| Species                     | Human                                                                                                                                                                                                                                                                                                                                                                                                              |
| Specificity/Sensitivity     | FLJ11151 - hypothetical protein FLJ11151,                                                                                                                                                                                                                                                                                                                                                                          |
| Immunogen                   | FLJ11151 (AAH06289, 1 a.a. - 314 a.a.) full-length human protein.<br>MSAAEAGGVFHRARGRTLDAFPAEKESEWKGPFFYFILGADPQFGLIKAWSTG<br>DCDNGGDEWEQEIRLTEQAVQAINKLNPKPKFFVLCGDLIHAMPGKPWRTEQT<br>EDLKRVLRVDRAIPLVLVSGNHDIGNTPTAETVEEFCRTWGDDYFSFWVGGV<br>LFLVLNSQFYENPSKCPSLKQAQDQWLDEQLSIARQRHCQHAIVFQHIPLFLESI<br>DEDDDYFFNLSKSTRKKLADKFIHAGVRVVFSGHYHRNAGGTYQNLDMMVSS<br>AIGCQLGRDPHGLRVVVVTAEKIVHRYYSLELSEKGIEDDLMDLIKKK |
| Notes                       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                                                 |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Applications                | Western Blot                                                                                                                                                                                                                                                                                                                                                                                                       |
| Recommended Dilutions       | Western Blot 1:500                                                                                                                                                                                                                                                                                                                                                                                                 |
| Application Notes           | Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.                                                                                                                                                                                                                                  |




## Images

Western Blot: CSTEP1 Antibody [H00055313-B01P] - Analysis of FLJ11151 expression in human liver.



Western Blot: CSTEP1 Antibody [H00055313-B01P] - Analysis of FLJ11151 expression in transfected 293T cell line by FLJ11151 polyclonal antibody. Lane 1: FLJ11151 transfected lysate(34.54 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00055313-B01P**

---

|                |                                                         |
|----------------|---------------------------------------------------------|
| HAF007         | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B        | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97019-5mg | Mouse IgG Isotype Control                               |
| NBP2-22723     | Recombinant Human CSTP1 His Protein                     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00055313-B01P](http://www.novusbio.com/reviews/submit/H00055313-B01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

