

Product Datasheet

PI4K2A 293T Cell Transient Overexpression Lysate H00055361-T02

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00055361-T02

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00055361-T02



H00055361-T02**PI4K2A 293T Cell Transient Overexpression Lysate**

Product Information	
Unit Size	0.1 ml
Concentration	mg/ml
Storage	Store at -80C. Avoid freeze-thaw cycles.
Buffer	1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue)

Product Description	
Description	Quality control test: Transient overexpression cell lysate was tested with Anti-PI4KII antibody by Western Blots. Plasmid: pCMV-PI4KII full-length
Gene ID	55361
Gene Symbol	PI4K2A
Species	Human
Specificity/Sensitivity	PI4KII 293T Cell Transient Overexpression Lysate(Denatured)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Lysate Type	293T Cell Transient Overexpression
Lysate Protein State	Denatured

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00055361-T02

H00055361-P01-2ug	Recombinant Human PI4K2A GST (N-Term) Protein
NBP2-13075	beta Amyloid Antibody (MOAB-2) - BSA Free
H00055361-D01P	PI4K2A Antibody
7268-CT-100	CTLA-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00055361-T02

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

