

Product Datasheet

SF20/MYDGF Antibody (2D4) - Azide and BSA Free H00056005-M03

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00056005-M03

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00056005-M03



H00056005-M03


SF20/MYDGF Antibody (2D4) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	2D4
Preservative	No Preservative
Isotype	IgG2a Kappa
Purity	IgG purified
Buffer	In 1x PBS, pH 7.4
Product Description	
Description	Novus Biologicals Mouse SF20/MYDGF Antibody (2D4) - Azide and BSA Free (H00056005-M03) is a monoclonal antibody validated for use in ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	56005
Gene Symbol	MYDGF
Species	Human
Immunogen	SF20/MYDGF (AAH10129, 1 a.a. ~ 173 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MAAPSGGWNGVGASLWAALLLGAVALRPAEAVSEPTTVAFDVRPGGVVHSFS HNVGPGDKYTCMFTYASQGGTNEQWQMSLGTSEDHQQHFTCTIWRPQGKSYL YFTQFKAEVARGAEIEYAMAYSKAAFERESDVPLKTEEFVTKTAVAHKPGAFKA ELSKLVIVAKASRTEL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	ELISA 1ng/mL, Immunocytochemistry/ Immunofluorescence 10ug/mL
Application Notes	Antibody reactive against recombinant protein on ELISA and Immunofluorescence. GST alone used as a negative control.




Images

Immunocytochemistry/Immunofluorescence: SF20/MYDGF Antibody (2D4) [H00056005-M03] - Analysis of monoclonal antibody to SF20/MYDGF on HeLa cell. Antibody concentration 10 ug/ml



ELISA: SF20/MYDGF Antibody (2D4) [H00056005-M03] - Detection limit for recombinant GST tagged SF20/MYDGF is approximately 1ng/ml as a capture antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00056005-M03

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP2-35172-100ug	Recombinant Mouse SF20/MYDGF Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00056005-M03

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

