

Product Datasheet

eIF4ENIF1 Antibody - Azide and BSA Free H00056478-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00056478-B01P

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00056478-B01P



H00056478-B01P

eIF4ENIF1 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description	
Host	Mouse
Gene ID	56478
Gene Symbol	EIF4ENIF1
Species	Human, Mouse
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 29395907).
Specificity/Sensitivity	Reacts with eukaryotic translation initiation factor 4E nuclear import factor 1.
Immunogen	EIF4ENIF1 (AAH33028.1, 1 a.a. - 985 a.a.) full-length human protein. MDRRSMGETESGDAFLDLKKPPASKCPHRYTKEELLDIKELPHSKQRPSCLSE KYDSGDGVDPEKWHASLYPASGRSSPVESLKKELDTDRPSLVRRIVDPRERV KEDDLDVVLSPQRRSFGGGCHVTAAVSSRRSGSPLEKDSGLRLLGGRRIGS GRIISARTFEKDHRLSDKDLRDLRDRERDFKDKRFRREFGDSKRVFGERRR NDSYTEEEPEWFSAGPTSQSETIELTGFDKILEEDHKGRKTRRRRTASVKEGI VECNGGVAEEDEVEVILAQEPAAADQEVPRDAVLPEQSPGDFDFNEFFNLDKVP CLASMIEDVLGEGSVSASRFSRWFSNPSRSGSRSSSLGSTPHEELERLAGLEQ AILSPGQNSGNYFAPIPLEDHAENKVDILEMLQKAKVDLKPILLSSLSANKEKLKE SSHSGVVLVSVEEVEAGLKGLKVDQVKNSTPFMAEHLEETLSAVTNNRQLKKD GDMTAFNKLVTMKASGTLPSQPKVSRNLESHLMSPAIEPGQPVPKNILQELL GQPVQRPASSNLLSGLMGSLIPTTSLGQRAPSPPLSQVFQTRAASADYLRPR IPSPIGFTPGPQQLLGDPFQGMKPMSPITAQMSQLELQQAALGLALPHDLAV QAANFYQPGFGKPQVDRTRDGFNRNRQQRVTKSPAPVHRGNSSSPAPAASITS MLSPSFTPTSVIRKMYESKEKSKEEPASGKAALGDSKEDTQKASEENLLSSSS VPSADRDSSPTTNSKLSALQRSSCSTPLSQANRYTKEQDYRPKATGRKTPTLA SPVPTTPFLRPVHQVPLVPHVPMVRPAHRLHPGLVQRMALAQGVHPQHLPSSL QTGVLPPGMDLSHLQGISGPILGQPFYPLPAASHPLLNPRTPLHLAMVQQQ LQRSVLHPPGSGSHAAVSVQTTQNVPSRSGLPHMHSQLEHRPSQRSSSPV GLAKWFGSDVLQQPLPSMPAKVISVDELEYRQ
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Dual RNAScope ISH-IHC
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Dual RNAScope ISH-IHC




Application Notes


This antibody is reactive against cell and transfected lysate in western blot, immunofluorescence, and as a detection antibody in ELISA. Dual ISH-IHC validation (PMID: 29395907).

Images


Western Blot: eIF4ENIF1 Antibody [H00056478-B01P] - Analysis of EIF4ENIF1 expression in Jurkat.




Immunocytochemistry/Immunofluorescence: eIF4ENIF1 Antibody [H00056478-B01P] - Analysis of purified antibody to EIF4ENIF1 on HeLa cell. (antibody concentration 10 ug/ml)



Western Blot: eIF4ENIF1 Antibody [H00056478-B01P] - Analysis of EIF4ENIF1 expression in Hela S3 NE.



Western Blot: eIF4ENIF1 Antibody [H00056478-B01P] - Analysis of EIF4ENIF1 expression in transfected 293T cell line by EIF4ENIF1 polyclonal antibody. Lane 1: EIF4ENIF1 transfected lysate(108.35 KDa). Lane 2: Non-transfected lysate.



Publications

Zahr SK, Yang G, Kazan H, A Translational Repression Complex in Developing Mammalian Neural Stem Cells that Regulates Neuronal Specification. *Neuron*. 2018-02-07 [PMID: 29395907] (Dual ISH-IHC, Mouse)

Yang G, Smibert CA, Kaplan DR, Miller FD. An eIF4E1/4E-T Complex Determines the Genesis of Neurons from Precursors by Translationally Repressing a Proneurogenic Transcription Program *Neuron*. 2014-11-06 [PMID: 25456498]

Amadei G, Zander MA, Yang G et al. A Smaug2-Based Translational Repression Complex Determines the Balance between Precursor Maintenance versus Differentiation during Mammalian Neurogenesis. *J Neurosci*. 2015-11-25 [PMID: 26609159]



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00056478-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
NBP1-89389PEP	eIF4ENIF1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00056478-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

