

# Product Datasheet

## eIF5A2 Antibody Pair H00056648-AP11

Unit Size: 1 Set

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00056648-AP11](http://www.novusbio.com/H00056648-AP11)

Updated 6/27/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00056648-AP11](http://www.novusbio.com/reviews/destination/H00056648-AP11)



**H00056648-AP11**

eIF5A2 Antibody Pair

| Product Information |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 1 Set                                                                                                                 |
| Concentration       | Concentration of individual antibodies may be found on the vial label. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                          |
| Clonality           | Monoclonal                                                                                                            |
| Preservative        | No Preservative                                                                                                       |
| Buffer              | 1x PBS, pH 7.4                                                                                                        |

| Product Description  |                                                                                                                                                                                                                                                                                                                                                                                      |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description          | <p>Quality control test: Standard curve using recombinant protein (H00056648-P01) as an analyte. This antibody pair set comes with matched antibody pair to detect and quantify protein level of human EIF5A2.</p> <p>Capture antibody: rabbit MaxPab affinity purified polyclonal anti-EIF5A2 (100 ug)<br/> Detection antibody: mouse monoclonal anti-EIF5A2 IgG1 Kappa (20 ug)</p> |
| Gene ID              | 56648                                                                                                                                                                                                                                                                                                                                                                                |
| Gene Symbol          | EIF5A2                                                                                                                                                                                                                                                                                                                                                                               |
| Species              | Human                                                                                                                                                                                                                                                                                                                                                                                |
| Notes                | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                   |
| Standard Curve Range | 3 ng/ml to 100 ng/ml                                                                                                                                                                                                                                                                                                                                                                 |
| Assay Type           | Sandwich ELISA                                                                                                                                                                                                                                                                                                                                                                       |

| Product Application Details |                                    |
|-----------------------------|------------------------------------|
| Applications                | Sandwich ELISA                     |
| Recommended Dilutions       | Sandwich ELISA 3 ng/ml - 100 ng/ml |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
[novus@novusbio.com](mailto:novus@novusbio.com)

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
[canada.inquires@bio-techne.com](mailto:canada.inquires@bio-techne.com)

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
[info.EMEA@bio-techne.com](mailto:info.EMEA@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [technical@novusbio.com](mailto:technical@novusbio.com)  
Orders: [orders@novusbio.com](mailto:orders@novusbio.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Antibody Pairs are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00056648-AP11](http://www.novusbio.com/reviews/submit/H00056648-AP11)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

