

Product Datasheet

TAPP1/PLEKHA1 293T Cell Transient Overexpression Lysate H00059338-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00059338-T01

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00059338-T01



H00059338-T01**TAPP1/PLEKHA1 293T Cell Transient Overexpression Lysate**

Product Information	
Unit Size	0.1 ml
Concentration	Please contact technical services for concentration.
Storage	Store at -80C. Avoid freeze-thaw cycles.
Buffer	1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue)

Product Description	
Description	Quality control test: Transient overexpression cell lysate was tested with Anti-PLEKHA1 antibody by Western Blots. Plasmid: pCMV-PLEKHA1 full-length
Gene ID	59338
Gene Symbol	PLEKHA1
Species	Human
Specificity/Sensitivity	TAPP1/PLEKHA1 293T Cell Transient Overexpression Lysate(Denatured)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Lysate Type	293T Cell Transient Overexpression
Lysate Protein State	Denatured

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00059338-T01

H00059338-Q01-10ug	Recombinant Human TAPP1/PLEKHA1 GST (N-Term) Protein
NB200-540	Complement C3 Antibody (11H9) - BSA Free
NBP1-86967	TAPP1/PLEKHA1 Antibody
AF887	AKT [p Ser473] Antibody - Pan Specific

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00059338-T01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

